

## Anti-GST Affinity Agarose

Cat. No: EA-IP-009

Size: 2 mL

**Note:** Do not centrifuge and use after mixing gently.

### Performance metrics

<b>Scope of application</b>	Affinity purification and immuno(co)precipitation of GST tag fusion proteins. The GST tag can be located at the N-terminus, C-terminus or in the middle of the protein, such as N-terminal GST fusion protein (GST - Protein), C-terminal GST fusion protein (Protein-GST) and Met-modified N-terminal GST fusion protein (Met-GST- Protein). Suitable for secreted proteins.
<b>Antibody properties</b>	Rabbit polyclonal antibody, IgG subtype 2a.
<b>Gel properties</b>	Agarose gel granules, average size 100~200 μm.
<b>Binding capacity</b>	1mL Sepharose 4B agarose particles, covalently conjugated to 6mg mouse-derived IgG. 1mL of affinity gel can purify or precipitate at least 1.2mg of GST fusion protein.
<b>Components</b>	1mL Anti-GST affinity gel in 1mL PBS with preservative and 50% glycerol.

### Matters Needing Attention

1. This product is only for scientific research by professionals and may not be used for clinical diagnosis or treatment.
2. For your safety and health, please wear a lab coat and disposable gloves.
3. This product is in the form of gel suspension, and the content of affinity gel is 50%. Gently re-suspend the gel suspension before use, and then use it as needed.
4. IP-WB samples are best prepared and used immediately to avoid affecting the experimental results.
5. Do not dry the gel, do not sonicate the gel, and do not allow acid treatment of the gel for more than 10 minutes.
6. The gel dosage in the usage method is a demonstration dosage prepared in a small amount. Please adjust the specific dosage according to the actual situation.

### Method of Application

#### 1. Sample Preparation of Target Proteins

##### 1) Sample processing serum and recombinant proteins

Collect serum or culture medium supernatant and detect the target protein concentration. If the target protein concentration is high, it is recommended to dilute it with 1×PBS to a final protein concentration of 10~100μg/mL for subsequent experiments.

##### 2) Sample processing of target protein for intracellular expression

- a. Blow off in case of adherent cells or take suspension cells from the cell culture flask and transfer them to a centrifuge tube, centrifuge at 1000 rpm for 5 min, and discard the supernatant.
- b. Re-suspend cells in 1× PBS pre-cooled at 4 °C, centrifuge at 1,000 rpm for 3 min, and discard the supernatant. Repeat once.
- c. Add the corresponding volume of cell lysate according to the amount of cells, and place on ice for 10~20 min after repeated pipetting.

Note: Generally, 1mL of cell lysis solution can process about 0.5~1×10<sup>7</sup> cells. To avoid degradation of your target protein, you can add protease inhibitors.

- d. Use a sonicator to treat the cell lysate until the cell lysate is transparent and no longer viscous. After placing on ice for 30 minutes, centrifuge at 12,000 rpm and 4°C for 10 minutes. Take the supernatant for subsequent experiments

### For Research Use Only

## 2. Column Installation and Incubation

### 1) Anti-GST Affinity Agarose preparation

- Gently re-suspend the Anti-GST Affinity Agarose, mix evenly, and take 40  $\mu$ L gel suspension (containing approximately 20  $\mu$ L gel) into a centrifuge tube.
- Add 10 times the gel volume (about 200  $\mu$ L) of 1 $\times$ PBS to gently resuspend the washed gel, centrifuge at 5 000 rpm for 30 seconds, discard the supernatant, and repeat this step three times.

**Note:** For multiple samples, the gel can be re-suspended and divided into several reaction tubes for separate reactions.

### 2) Binding of target protein to Anti-GST Affinity Agarose

- Incubation: Add 200  $\mu$ L of the prepared sample to the washed gel, and incubate on a shaker at room temperature for 2 hours. It can also be incubated at 4°C overnight or longer.
- Washing: After incubation, centrifuge at 5000rpm for 30 seconds and discard the supernatant. Add 200  $\mu$ L 1 $\times$ PBST, mix gently, wash the gel, centrifuge at 5000 rpm for 30 seconds, discard the supernatant, and repeat this step 4 times.

### 3) Target protein elution

This instruction manual provides the following two target protein elution schemes. Please choose different target protein elution methods according to the needs of later detection.

#### Denaturing elution method

This method is only suitable for SDS-PAGE detection.

- Add 16  $\mu$ L 1 $\times$ PBS and 4  $\mu$ L 5 $\times$  loading buffer, boil the sample for 5 minutes, cool it down to room temperature and centrifuge.
- Take the supernatant and run SDS-PAGE in preparation for subsequent Western Blot detection.

#### Acid elution method

Acidic elution method has low cost, short operational time, generally does not cause protein denaturation, and facilitates subsequent analysis and detection of proteins.

- Add pre-cooled acidic eluent pH 3.0, 10 times of the gel volume (approximately 200  $\mu$ L), to the above precipitate, suspend the affinity gel, and incubate at room temperature for 5 minutes.  
**Note:** An acidic environment will shorten the service life of the gel. The contact time between the gel and the acidic eluent should be shortened as much as possible. It is recommended not to exceed 10 minutes.
- After the incubation, centrifuge at 5000 rpm for 30 seconds at 4°C, transfer the supernatant to a new centrifuge tube, and immediately add 1/10 volume of neutralizing solution pH 8.0 and mix well. The supernatant is the eluted GST-tagged protein.
- Process and store proteins according to subsequent experimental needs.

#### Gel washing and regeneration

If the Affinity Agarose needs to be reused, it must be washed and regenerated immediately after elution.

- Wash once with 10 times the gel volume of acidic eluent, 10 times the gel volume of neutralizing solution, and 10 times the gel volume of 1 $\times$ PBS.
- Wash once more with 3 times the volume of PBS containing preservative and 50% glycerin.
- Store in an equal volume of gel containing PBS, preservatives and 50% glycerol, and store sealed at -20°C.

## For Research Use Only

## Background

Anti-GST Affinity Agarose is made by covalently conjugating high-quality GST tag antibodies to agarose gel. It has the characteristics of high loading capacity, high specificity, stable properties, and can be used repeatedly, and can be used for affinity purification, immunoprecipitation (IP), co-immunoprecipitation (Co-IP) and other immunoprecipitation related experiments.

## Storage

-20°C for 12 months.