



The professional cell culture empowers a healthler work

Nanobacteria Removal Complete Medium (10) (DMEM (High glucose))

Cat. No.: PM150210B-HR

Size: 125mL×4

General Information

Product Form Liquid

Concentration Ready-to-use

Components DMEM (High glucose)[PM150210]+10% Nutrient+1% Supplement+Anti-

Nanobacteria Treatment Reagent[P-CMR-002]

Bacterial detectionNegativeFungal detectionNegativeMycoplasmal detectionNegativeEndotoxin level< 3 EU/mL</th>ShippingIce bag

Storage 2-8°C, Shading Light

Expiry date 6 months

Product Introduction











DMEM (Dulbecco's Modified Eagle Medium) was developed on the basis of MEM medium. Compared with MEM medium, the content of amino acid increased by 2 times, the content of vitamin increased by 4 times, and the content of non-essential amino acid, trace iron ion and sodium pyruvate were increased by 4 times.

The glucose content of DMEM medium was originally designed as 1000 mg/L (low Glucose type), and then developed into 4500 mg/L (high Glucose type), which has been widely used in cell culture.

DMEM (High glucose) was widely used in fast growth, low adhesion cells, hybridoma myeloma cells, clone cells, DNA transfected transformation cells, various primary virus host cells, single cell culture and vaccine production.

DMEM (High glucose) contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors, so the product should be used with serum or serum-free additives.

Nanobacteria and their decomposition complexes are the common contaminant in cell cultures that co-exists with cells. Antibiotics are usually ineffective. Nanobacteria grows competitively with cells, which is unfavorable to cell growth, and in severe cases causes cell death. At present, many cells are contaminated with nanobacteria, which has a great impact on cell culture and subsequent experiments. The common characteristics of cells contaminated by nanobacteria are as follows: (1) The medium is not turbid, but when the cells are observed under a microscope, there are many "small black spots" around the cells or in the culture medium. With the extension of culture time, the "small black spots" gradually increase, and they cannot be removed by changing the culture medium or washing the cells.

- (2) The cells contaminated by the "small black spots" consume fast nutrients and require frequent replacement of the culture medium.
- (3) Nanobacteria-contaminated cells grow slowly, have poor cell states, and are severely vacuolated. They may even cause changes in cell morphology. Therefore, it is of great significance to clean up nanobacteria contamination in cell culture.

Anti-Nanobacteria media is a new generation product developed by our team on the basis of Biomocin products. It contains special ingredients to remove nanobacteria. This product has been tested on hundreds of cells and verified by long-term experiments. It is harmless to cells and has a significant effect on removing and inhibiting nanobacteria.

Guidelines for use



- 1. Pricella's cell culture media undergoes strictquality control to ensure sterility, but may get contaminated during use. Follow these guidelines for sterile handling toavoid contamination.
- 2. Always wipe your gloved hands and work area with 70% ethanol.
- 3. Wipe the outside of the containers, flasks, plates, and dishes with 70% ethanol before placing them in the cell culture hood.
- 4. Use sterile pipette tips and pipettes to work with liquids, and use each pipette tip only once to avoid cross contamination. Do not unwrap sterile pipettes until they are ready to beused. Keep pipettes and tips within the clean work area.
- 5. Do not talk while performing sterile procedures and perform your cell culture as efficiently and carefully as by Elabscier possible tominimize contamination.

Quality control

Standard evaluations for cell culture media are pH, osmolality, endotoxins and sterility testing for liquid products, cell growth experiments.


Notes

Tel:1-832-243-6086 Fax:1-832-243-6017 Toll-free:1-888-852-8623 Email: techsupport@pri-cella.com Web: www.pri-cella.com





- 1. This product is for research use only.
- 2. This product is sterilized by 0.1 μm filtration.
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.
- 4. It is not suitable for long time storage at room temperature.
- 5. This product is a ready-to-use medium. If there is no special need, don't add serum, penicillin and streptomycin. It can be used directly.









Tel:1-832-243-6086 Fax:1-832-243-6017 Toll-free:1-888-852-8623 Web: www.pri-cella.com