EasySort™ Human Naïve Pan T Cell Isolation Kit

Cat. No: MIH006N Size: 10/100/200 Assays

Component	Component Name				Storage
MIH006NA	EasySort™ Human Naïve Pan T Beads Streptavidin 1.0-N	300 µL	1.0 mL×3	1.0 mL×6	2-8°C
MIH006NB	EasySort™ Human Naïve Pan T Cell Isolation Cocktail	375 μL	1.25 mL×3	1.25 mL×6	2-8°C
Manual		1 сору			

Storage

Store at 2-8°C with shading light for 1 year. Avoid freezing and thawing.

Description

Human Naïve Pan T cell isolation kit utilizes a negative selection strategy to isolate naïve pan T cells from fresh or frozen human PBMC sample. The principle of this kit is to use different biotinylated monoclonal antibodies labeling non-Naïve Pan T cells, followed by streptavidin conjugated magnetic beads incubation. Unwanted labeled cells are efficiently removed by magnetic isolation, and high purity Naïve Pan T cells are isolated.

EasySort™ Human Naïve Pan T Cell Isolation Kit can help researcher isolate high purity human Naïve Pan T cells with simple experimental procedure. The kit is suitable for isolation of Naïve Pan T cells from fresh human PBMC or frozen PBMC, and the isolated Naïve Pan T cells can be directly used for downstream applications. The Naïve Pan T cells isolated from normal PBMC using this kit is typically 94.2 ± 1.6%.

Reagents and Materials Not Supplied

1. Reagents:

PBS, fetal bovine serum (FBS), EDTA, Human peripheral blood mononuclear cells separation solution, DNase I

2. Materials:

70 µm mesh nylon strainer, 1.5 mL/2 mL EP tube, 15 mL/50 mL centrifuge tube, flow tube

3. Instrument:

Optical microscope, centrifuge, 5 mL magnetic rack

For Research Use Only



Experimental Operation

NOTE: The following operations must be performed under sterile conditions

Isolation buffer preparation

Add fetal bovine serum (final concentration of 2%) and EDTA (final concentration of 2 mM) to PBS buffer and filter the prepared buffer with 0.22 µm filter.

NOTE: Sealed store the prepared buffer at 4°C and use within 1 week. In addition, 2% fetal bovine serum can be replaced by 0.5% BSA.

> Sample Preparation and Processing

1. Fresh human PBMC: PBMC sample is obtained from fresh human whole blood by density gradient centrifugation. Wash PBMC twice with isolation buffer, centrifuge at 300 g for 5 min, filter the PBMC through a 70 μ m mesh nylon strainerand and adjust the cell density to 1×10^8 cells/mL for cell isolation.

Note: The best separation effect can be achieved when the freshly collected human blood is separated within 1 hour. Approximately 1×10⁷ PBMC can be obtained from 10 mL of human blood.

2. Frozen PBMC: incubate the frozen PBMC should be incubated with DNase I solution (PBS) at a concentration of 100 μg/mL for 15 min at room temperature before cell isolation. Wash sample twice with isolation buffer, centrifuged at 300 g for 5 min. Filter aggregated suspensions through a 70 μm mesh nylon strainer and adjust cell density at 1×10⁸ cells/mL.

Cell Isolation

a) Prepare 100 μ L of cell suspension (about 1×10⁷ cells), add 37.5 μ L Human Naïve Pan T Cell Isolation Cocktail, mix fully and incubate for 5 min at room temperature.

Note: Please make sure the cells are single-cell suspension.

- b) Add isolation buffer to a final volume of 2 mL, centrifuge at 300 g for 5 min. Discard the supernatant, and then resuspend the cells with 100 µL isolation buffer.
- c) Wash Beads Streptavidin 1.0-N: Vortex beads for 20 seconds, add 30 µL Beads in 1.5 mL EP tube. Put the tube on a 5 mL magnetic rack (self-provided) and stand for 30 seconds. Remove the supernatant, then resuspend beads with 1 mL isolation buffer, and stand for 5 minutes at room temperature. Remove the supernatant, then resuspend beads with 30 µL isolation buffer.
- d) Transfer the cells to the bottom of the flow tube (Note: Avoid adding along tube walls), add 30 μL washed Human Naïve Pan T Beads Streptavidin 1.0-N, mix gently and incubate at room temperature for 5 min.


Note:

If more than 1×10⁷ cells are to be isolated, increase the amount of Human Naïve Pan T Cell Isolation Cocktail and Human Naïve Pan T Beads Streptavidin 1.0-N proportionally while ensuring the cell density remains 1×10⁸ cells/mL. If fewer than 1×10⁷ cells are to be isolated, resuspend the cells with 100 μL isolation buffer, add 37.5 μL Human Naïve Pan T Cell Isolation Cocktail and 30 μL washed Human Naïve Pan T Beads Streptavidin 1.0-N.

For Research Use Only

- ♦ The 5 mL flow tube is suitable for less than 1×10⁸ cells.
- e) Add isolation buffer to a final volume of 2.5 mL, mix gently with a pipette by blowing up and down for 7-8 times until no particles of magnetic beads are visible. Put the tube on a 5 mL magnetic rack (self-provided) and stand for 5 min.
 - Note: Please mix the liquid thoroughly to avoid the magnetic beads clumping and affecting the isolation efficiency.
- f) Transfer the cell suspension to a clean centrifuge tube, centrifuge at 300 g for 5 min. Discard the supernatant, resuspend the cells with buffer required for the subsequent experiments.

Typical data



In the above example, the purities of human naïve pan T cells (CD3+CD45RA+CD45RO-CD197+) in normal human peripheral blood mononuclear cell samples was 6.0% before sorting and 95.2% after sorting.





Fluorochrome-conjugated antibody	Cat.	Company	
Elab Fluor® Violet 450 Anti-Human CD3	E AD E10010	Elabscience	
Antibody[OKT-3]	E-AB-F1001Q		
PE Anti-Human CD45RA Antibody[HI100]	E-AB-F1052D	Elabscience	
APC Anti-Human CD45RO Antibody[UCHL1]	E-AB-F1139E	Elabscience	
FITC Anti-Human CD197/CCR7 Antibody[G043H7]	E-AB-F1159C	Elabscience	

Cautions

- 1. This kit is for research use only.
- 2. Please take safety precautions and follow the procedures of laboratory reagent operation.
- 3. Avoid freezing and thawing during the use and storage of the beads.
- 4. The data presented for this product are based on testing of peripheral blood mononuclear cells (PBMCs) from healthy individuals. When processing peripheral blood from non-healthy individuals, to ensure optimal experimental outcomes, first assess the naïve pan T-cell content in the sample. If the naïve pan T-cell content is below 6%, it is recommended to perform a pre-sorting test, appropriately increase the amount of antibodies and beads per assay, or conduct a secondary sorting step.
- 5. The cell clusters in the cell suspension will affect the purity of cell isolation. Therefore, cell suspension should be filtered with a 70 µm mesh nylon sieve before formal isolation.
- 6. Cell suspension should be isolated immediately after preparation, the longer the storage time, thegreaterthe impact on cell activity.
- 7. The cell suspension and reagents should be added directly to the bottom of flow tube to avoid sticking to the wall, resulting in insufficient reaction and affecting the isolation efficiency.
- 8. In order to ensure the activity of the cells, the whole process of the experiment should be completed on ice as much as possible, except for the incubation at room temperature.
- 9. It is recommended to use low adsorption pipette tips and centrifuge tubes to avoid the loss of magnetic beads and antibodies due to adsorption.
- 10. The kit should be used in combination with a magnetic rack.

For Research Use Only