

Anti-GST immunomagnetic beads

Cat. No: EA-IP-009M

Size: 1 mL

Note: Do not centrifuge and use after mixing gently.

Performance metrics

Scope of application	Applied to the immunoprecipitation of fusion proteins or protein complexes with GST tags. The GST tag can be located at the N-terminus, C-terminus or in the middle of the protein, such as N-terminal GST fusion protein (GST - Protein), C-terminal GST fusion protein (Protein-GST) and Met-modified N-terminal GST fusion protein (Met-GST- Protein).
Antibody properties	GST-Tag Mouse Ab: Mouse IgG.
Gel properties	Agarose gel granules, average size 100~200 μm.
Binding capacity	1 mL of magnetic bead suspension containing 20 mg of magnetic beads covalently conjugated to approximately 1 mg of Anti-GST mouse monoclonal antibody. 1mL of Anti-GST immunomagnetic beads can precipitate 1~2 mg of GST fusion protein.
Components	0.25mL Anti-GST immunomagnetic beads, stored in 0.75mL PBS containing preservatives.

Matters Needing Attention

1. This product is only for scientific research by professionals and may not be used for clinical diagnosis or treatment.
2. For your safety and health, please wear a lab coat and disposable gloves.
3. This product provides affinity magnetic beads in the form of suspension. Gently re-suspend the magnetic bead suspension before use, and then use it as needed.
4. Do not centrifuge, freeze or dry the magnetic beads, do not sonicate the magnetic beads, and do not allow acid treatment of the magnetic beads for more than 10 minutes.
5. When mixing the magnetic beads, please use methods such as gentle pipetting with a pipette, gentle vortexing, inversion, and shaker mixing. Do not use sonication and other methods.
6. The relevant reagents used must be prepared by the laboratory.

Method of Application

1. Sample Preparation of Target Proteins

1) Sample processing serum and recombinant proteins

Collect serum or culture medium supernatant and detect the target protein concentration. If the target protein concentration is high, it is recommended to dilute it with 1×PBS to a final protein concentration of 10~100μg/mL for subsequent experiments.

2) Sample processing of target protein for intracellular expression

- a. Blow off in case of adherent cells or take suspension cells from the cell culture flask and transfer them to a centrifuge tube, centrifuge at 1000 rpm for 5 min, and discard the supernatant.
- b. Re-suspend cells in 1× PBS pre-cooled at 4 °C, centrifuge at 1,000 rpm for 3 min, and discard the supernatant. Repeat once.
- c. Add the corresponding volume of cell lysate according to the amount of cells, and place on ice for 10~20 min after repeated pipetting

Note: Generally, 1mL of cell lysate can process about 0.5~1×10⁷ cells. To avoid degradation of the target protein, you can add a protease inhibitor.

For Research Use Only

- d. Use a sonicator to treat the cell lysate until the cell lysate is transparent and no longer viscous. After placing on ice for 30 minutes, centrifuge at 12,000 rpm and 4°C for 10 minutes. Take the supernatant for subsequent experiments.

2. Column Installation and Incubation

1) Anti-GST Immunomagnetic beads preparation

- a. Gently re-suspend the Anti-GST immunomagnetic beads, mix evenly, and take 40 µL of the magnetic bead suspension (containing approximately 10 µL of magnetic beads) into a centrifuge tube.
- b. Add 500 µL of 1×PBS to gently re-suspend and wash magnetic beads, let stand on the magnetic stand for 10 seconds, discard the supernatant, and repeat the above steps twice.

Note: For multiple samples, the magnetic beads can be re-suspended and divided into several reaction tubes for separate reactions.

2) Binding of target protein to Anti-GST immunomagnetic beads

- a. Incubation: Add 500 µL of the prepared sample to the washed magnetic beads, and incubate on a shaker at room temperature for 2 hours. It can also be incubated at 4°C overnight or longer.
- b. Washing: After incubation, perform magnetic separation and discard the supernatant. Add 500 µL 1×PBST, mix gently, wash the magnetic beads, magnetically separate, and discard the supernatant. Repeat 3 times.

3) Target protein elution

This instruction manual provides the following two target protein elution schemes. Please choose different target protein elution methods according to the needs of later detection.

Denaturing elution method

This method is only suitable for SDS-PAGE detection.

- a. Add 20 µL 1×PBS and 5 µL 5× loading buffer, boil the sample for 5 minutes, cool it down room temperature and centrifuge.
- b. Take the supernatant and run the SDS-PAGE in preparation for subsequent Western Blot detection.

Acid elution method

Acidic elution method has low cost, short operation time, generally does not cause protein denaturation, and facilitates subsequent analysis and detection of proteins.

- a. Add pre-cooled acid eluent pH 3.0, 0.5 mL or 20 times the volume of magnetic beads, to the above precipitation, suspend the magnetic beads, and incubate at room temperature for 5 minutes.

Note: Acidic environment will shorten the service life of immunomagnetic beads. The contact time between magnetic beads and acidic eluent should be shortened as much as possible. It is recommended not to exceed 10 minutes.

- b. After the incubation, magnetically separate, transfer the supernatant to a new centrifuge tube, and immediately add 1/10 volume of pH 8.0 neutralizing solution and mix well.
- c. Process and store proteins according to subsequent experimental needs.

Background

Anti-GST immunomagnetic beads are made of high-quality GST antibodies covalently conjugated to magnetic beads. They can specifically bind to GST labeled proteins in the sample, thus being used for immunoprecipitation or co precipitation of GST labeled fusion proteins or their protein complexes.

Storage

For Research Use Only

-20°C for 12 months.