

## OrgaMatrigel Organoid Culture Specialized Matrix

Cat. No. : PB180536

Size: 5mL/5mL×2

### General Information

<b>Product Form</b>	Liquid state after thawing at 4°C
<b>Concentration</b>	Protein concentration: 8–12 mg/mL
<b>Size</b>	5mL/5mL×2
<b>Matrix Components</b>	Laminin, Type IV collagen, heparan sulfate proteoglycans (HSPG), Entactin, and various growth factors
<b>Bacterial detection</b>	Negative
<b>Fungal detection</b>	Negative
<b>Mycoplasma detection</b>	Negative
<b>Organoid Growth Assay</b>	Various organoids exhibit healthy growth and normal morphology
<b>Endotoxin level</b>	< 4 EU/mL
<b>Storage</b>	≤ -15°C
<b>Shipping</b>	Ice bag
<b>Expiration date</b>	24 months

### Background

OrgaMatrigel Organoid Culture Specialized Matrix is a natural basement membrane matrix extracted from mouse tumors rich in extracellular matrix proteins, containing laminin (Laminin), type IV collagen (Col-IV), heparan sulfate proteoglycans (HSPG), entactin (Entactin), growth factors, and other components.

This product has been meticulously optimized by Pricella's Technology team as a high-performance matrix gel specifically designed for in vitro organoid culture. It simulates the biochemical and physical characteristics of the in vivo extracellular matrix (ECM), providing an ideal 3D growth microenvironment for stem/progenitor cells, tumor cells, or other primary cells. Through extensive testing, this product has demonstrated its ability to support efficient organoid formation, expansion, and functional maintenance.

### Note

1. Aseptic technique must be maintained during product use to avoid contamination.
2. Avoid storing the matrix gel in a frost-free refrigerator. If the matrix gel will be used within 1 year and requires frequent access, thaw and aliquot it for storage to minimize repeated freeze-thaw cycles. Store at -20°C under normal conditions. For long-term storage (beyond 1 year) or experiments requiring extremely high stability, store at -80°C.
3. The matrix gel solidifies into a gel state at 37°C and will begin to solidify when the temperature is higher than 10°C. To ensure optimal performance, pre-cool pipette tips and perform all operations on ice.
4. All products must be used within their shelf life. Discard immediately if expired.