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1. Tseng S.Ja. An acid degradable, lactate oxidizing nanoparticle formulation for non-small cell lung cancer virotherapy[J]. *Nano Today*. IF:18.962
2. Salman T M, Iyanda M A, Alli-Oluwafuyi A M, et al. Telfairia occidentalis stimulates hepatic glycolysis and pyruvate production via insulin-dependent and insulin-independent mechanisms[J]. *Metabolism Open*, 2021, 10(1-10):100092. IF:8.694
3. Wang H, Huang Q, Zhang Z, et al. Transient post-operative overexpression of CXCR2 on monocytes of traumatic brain injury patients drives monocyte chemotaxis toward cerebrospinal fluid and enhances monocyte-mediated immunogenic cell death of neurons in vitro[J]. *Journal of Neuroinflammation*, 2022. IF:7.573
4. Carafa V, Russo R, Della Torre L, et al. The Pan-Sirtuin Inhibitor MC2494 Regulates Mitochondrial Function in a Leukemia Cell Line[J]. *Frontiers in Oncology*, 2020, 10: 820. IF:4.848
5. Sun J, Leng P, Li X, et al. Salvianolic acid A promotes mitochondrial biogenesis and mitochondrial function in 3T3-L1 adipocytes through regulation of the AMPK-PGC1 α signalling pathway. *Adipocyte*. 2022;11 (1):562-571. IF:3.553
6. Xu X, Cui Y, Li C, et al. SETD3 Downregulation Mediates PTEN Upregulation-Induced Ischemic Neuronal Death Through Suppression of Actin Polymerization and Mitochondrial Function. *Mol Neurobiol*. 2021; 58 (10):4906-4920. IF:3.464
7. Gao S, Li N, Chen R, et al. Bushen Huoxue Formula Modulates Autophagic Flux and Inhibits Apoptosis to Protect Nucleus Pulposus Cells by Restoring the AMPK/SIRT1 Pathway. *Biomed Res Int*. 2022; 2022:8929448. IF:3.047
8. Cheng F, Yu J, Zhang X, et al. CircSEC31A Promotes the Malignant Progression of Non-Small Cell Lung Cancer Through Regulating SEC31A Expression via Sponging miR-376a[J]. *Cancer Management and Research*, 2020, Volume 12:11527-11539. IF:2.886
9. Soebagjo H D, Fatmariyanti S, Paulus Sugianto N, et al. Detection of the Calcium and ATP Role in Apoptosis of Retinoblastoma Culture Cells through Caspase-3 Expression[J]. *Research Journal of Pharmacy and Technology*, 2019, 12(03): 1307-1314.
10. Al-Ziaydi A G, Al-Shammari A M, Hamzah M I, et al. Newcastle disease virus suppress glycolysis pathway and induce breast cancer cells death[J]. *VirusDisease*, 2020: 1-8.