

One-component TMB Substrate

Catalog No: E-IR-R201

30 mL/100 mL/500 mL/1000 mL /2000 mL

This manual must be read attentively and completely before using this product. If you have any problems, please contact our Technical Service

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Email: techsupport@elabscience.com

Website: www.elabscience.com

Please kindly provide us the lot number (on the outside of the box) of the kit for more efficient service.

Description

TMB is a popular chromogenic substrate for HRP detection in ELISA and is available in several formats. The soluble One-component TMB Substrate are supplied as Ready-to-Use mono-component liquids that require no preparation before use. The product mainly contains 3,3',5,5'-tetramethylbenzidine and contains no DMF or DMSO.

One-component TMB Substrate yields a blue color ($A_{\text{max}} = 370 \text{ nm}$ and 652 nm) when reacted with peroxidase in microwell applications. The product changes to yellow ($A_{\text{max}} = 450 \text{ nm}$) upon addition of a sulfuric or phosphoric acid stop solution.

Application

Peroxidase-based enzyme immunoassays ELISA is recommended.

Storage

This product is stable for 1 years stored at $2^{\circ}\text{C} \sim 8^{\circ}\text{C}$. Shading light.

Product features

- > **Economical:** Cost less compared with similar products.
- > **Convenience:** mono-component, ready-to-use.
- > **High sensitivity:** Save the dosage of antibody.
- > **Strong stability:** Valid for 1 years at $2 \sim 8^{\circ}\text{C}$.
- > **Stable signal:** The signal is continuous steady after adding stop solution.
- > **Safe:** Noncarcinogenic.

Assay procedure

1. Add $100 \mu\text{L}$ of the TMB Substrate Solution to each microplate well.
2. Incubate for $5 \sim 30 \text{ min}$ according to the reaction system.
3. Stop reaction by adding $50 \mu\text{L}$ stop solution (1M HCl or $1\text{M H}_2\text{SO}_4$).
4. Measure the absorbance of each well at 450 nm .

Notes

1. This substrate is light sensitive, avoid direct sunlight during operation.
2. This substrate is sensitive to contamination of oxidizing agent. Avoid oxidizing agent contamination during operation.
3. Never pipette directly from the bottle. Pour out required amount into a tube and pipette from the tube.
4. Do not return excess TMB to primary storage container.