

## Aptplex™ Human Th1/Th17 12-Plex Panel

Catalog No: MPA0011 ~~Aptplex™ Human COVID-19~~

### 8-Plex Panel

Catalog No: MPA001

Product size: 96 T

### Intended Use

Aptplex™ Human Th1/Th17 12-Plex Panel ~~Human COVID-19~~

~~8-plex Panel~~ is based on multiplex bead-based technology,

enabling simultaneous quantification of multiple analytes from a

single sample. This kit is suitable for the *in vitro* quantitative

detection of concentrations of the following cytokines in human

serum, plasma, cell culture supernatants, and other biological

fluids:

IL-17A, IL-6, TNF- $\alpha$ , IFN- $\gamma$ , IL-12p70, IL-10,

CCL5(RANTES), CXCL10(IP-10), IL-1 $\beta$ , IL-18,

sCD54(ICAM-1), IL-2RA(sCD25)IL-6, IL-1 $\beta$ , TNF- $\alpha$ , IFN- $\gamma$ ,

CXCL10 (IP-10), CCL2 (MCP-1), IL-8 (CXCL8), and IL-10.

### Components

Component	Component Name	96T	Storage
MPA0011A	Premixed Antibody-Conjugated Beads	2.4 mL×2	2-8°C <del>shading</del> Protected from light
MPA0011B	Biotinylated Detection Antibodies	4.8 mL×2	2-8°C
MPA0011C	SA-PE (ready to use)	4.8 mL×2	2-8°C <del>shading</del> Protected from light
MPA0011D	Lyophilized Standard	2 vials	2-8°C
MPA0110E	Assay Buffer	5 mL×1	2-8°C
MPA0011F	Wash Buffer	30 mL×2	2-8°C
	Plate Sealing Film	5 pieces	

### Introduction

Aptplex™ Human Covid-19 8-Plex Panel is based on multiplex

bead-based technology, enabling simultaneous quantification of

multiple analytes from a single sample. This kit is suitable for the

*in vitro* quantitative detection of concentrations of IL-6; IL-1 $\beta$ ;

TNF- $\alpha$ ; IFN- $\gamma$ ; CXCL10(IP-10); CCL2(MCP-1); IL-

8(CXCL8); IL-10 cytokines in human serum, plasma, cell

culture supernatants, and other biological fluids.

### Detection Principle

The Aptplex™ assay is a multiplex bead based immunoassay that

utilizes ~~uses~~ antibody-conjugated magnetic beads with distinct

~~mean~~ fluorescence intensities to capture target antigens from

~~samples~~ simultaneously. When selected capture beads are mixed

and incubated with samples containing specific ~~Each~~ targets

antigen is recognized by the ~~a specific~~ capture antibodies, each

analyte binds to its corresponding on the bead and ,

Subsequently ~~binds with a corresponding~~ biotinylated detection

antibody ~~to form a ies~~ are added, each of which binds to a specific

analyte, forming a "sandwich" complex of bead-analyte-detection

antibody "sandwich complex". Streptavidin-phycoerythrin (SA-

PE) ~~is then introduced to binds to~~ the biotinylated detection

antibodies producing a fluorescent signal proportional to ~~The~~

~~fluorescence intensity specific to the amount of~~ each analyte. The

~~fluorescence of each bead~~ is measured using flow cytometry. By

and correlating the results with a standard curve generated from

antigen standards, quantitative detection of target proteins in

~~samples is achieved to determine analyte concentrations.~~

### Detection Sample Types

☒ sSerum ☒ EDTA Pplasma

☒ eCell culture supernatants ☒ oOther biological fluids

### Storage

Store protected from light at 2-8 °C with a shelf life of 12 months.

After opening, protect from light and store at 2-8 °C, stable for up

to 30 days. Reconstituted Standard should be used within 24

hours.

Material	Storage Conditions	Stability / Notes
Unopened kit	2-8 °C, protected from light	12 months
Opened kit	2-8 °C, protected from light	Up to 30 days
Reconstituted standard	2-8 °C, protected from light	Use within 24 hours

### Sample Collection and Materials Not Supplied

#### 1) Serum

Allow ~~samples~~ whole blood to clot for 1 hour at room temperature

or overnight at 2-8 °C, ~~before centrifugation then centrifuge~~ for 20

min at 1000 ×g at 2-8 °C. Collect the supernatant ~~to carry out for~~

the assay.

#### 2) Plasma

Collect ~~plasma~~ using EDTA-Na<sub>2</sub> as an anticoagulant. Centrifuge

~~samples for~~ 15 min at 1000 ×g at 2-8 °C within 30 min of

collection. Collect the supernatant ~~to carry out for~~ the assay.

#### 3) Cell culture supernatant or other biological fluids

批注 [S1]: Please confirm the correct trademarked product name and its capitalization, should it be 'Aptplex', 'APTPlex', or 'APTplex'? Consistent use of the correct capitalization is important not only for clarity and professionalism in the manual but also to accurately reflect the official trademark as registered by the company. And once you confirm it, revise it throughout the document

带格式的: 左, 段落间距段后: 0.5 行

带格式的: 字体: 倾斜

带格式的: 字体: 六号

带格式的: 字体: 六号

带格式的: 字体: 六号

批注 [S3]: A table format makes storage conditions clear, concise, and easy to reference.

带格式的: 字体: 六号

带格式的: 字体: 六号

带格式的: 下标

批注 [S2]: The current Detection Principle reads more like a step-by-step protocol. It should be revised to focus on the scientific mechanism of the assay rather than procedural details, making it concise and clear for the manual. I strongly suggest to keep my version.

Centrifuge samples for 20 min at 1000 ×g at 2-8 °C. Collect the supernatant to carry out for the assay.

#### 4) Instrument Materials Not Supplied

- U-bottom 96-well transparent plates;
- Vortex mixer;
- Incubators suitable for 96-well plate;
- Magnetic separator;
- Flow cytometry (with PE , APC and APC/Cy7 detection channels).

#### Standard Preparation Procedure

1. Prepare eight 0.6 mL microcentrifuge tubes as standard dilution tubes and label them 0-7. Leave tube 7 empty. Add 150 µL of Assay Buffer to tubes 0-6.

#### 2. Reconstitute the lyophilized standard

- The standard is supplied as a lyophilized powder. Before reconstitution, briefly centrifuge the vial tube at 500 ×g for 10 seconds to collect the powder at the bottom of the tube.
- Add 500 µL of Assay Buffer to the vial. Let it stand for 5 min,
- then gently pipette up and down 2-3 times to ensure the Mix gently until the lyophilized standard is completely dissolved, and mixed.
- Transfer the entire solution to tube 7. This is the highest concentration standard solution.

#### 3. Prepare serial dilutions

- Pipette/Transfer 50 µL from tube 7 into tube 6 and mix thoroughly (1:3.2 dilution). This is the 1:3.2 dilution standard solution.
- Serially dilute the standard solution by transferring 100 µL from each tubes 6 to 2 in to the next respective lower tube (tubes 5 to 1: 6 → 5 → 4 → 3 → 2 → 1), mixing thoroughly at each step. This creates the standard solutions for tubes 5-1.
- The Assay Buffer in Tube 0 contains Assay Buffer only and serves as the zero concentration standard.

Note: The specific exact concentration of the highest concentration standard solution may vary. Please refer to is indicated on the standard vial label for details.

#### Procedure

##### 1. Bead and Sample incubation

1. Add 50 µL of Premixed Antibody-Conjugated Beads (vortex for at least ≥ 15 seconds before adding use) and 50 µL of sample or standard to each well of a 96-well plate. Cover the plate with a sealing film and incubate on a microplate shaker at room temperature, 600 rpm, protected from light, at 600 rpm for 1 hour.

##### 2. Detection antibody incubation

2. After incubation, place the plate on a magnetic separator for 1 minute. Decant and remove the supernatant, then add 100 µL of Biotinylated Detection Antibodies to each well. Cover the plate, with a sealing film and incubate under same conditions on a microplate shaker at room temperature,

protected from light, at 600 rpm for 1 hour.

#### 3. SA-PE incubation and washing

3. After incubation, place the plate on a magnetic separator for 1 minute. Decant/remove the supernatant, and wash once with add 200 µL of Wash Buffer to each well, place the plate on the magnetic separator for 1 minute, and decant the supernatant. Remove the wash buffer. Then, add 100 µL of SA-PE to each well. Cover the plate with a sealing film and incubate on a microplate shaker at room temperature, 600 RPM, protected from light, at 600 rpm for 30 minutes. 0.5 hours. Following incubation, wash twice with 200 µL of Wash Buffer, separating each time.

#### 4. Final resuspension and detection

- After incubation, place the plate on a magnetic separator for 1 minute. Decant the supernatant, add 200 µL of Wash Buffer to each well, place the plate on the magnetic separator for 1 minute, and decant the supernatant. Finally, add Resuspend the beads in 200 µL of Wash Buffer to each well to resuspend the beads.

and Proceed with detection using the a flow cytometer under appropriate instrument settings.

(The experimental protocol conditions may be adjusted by the user according to their specific experimental conditions and sample type or available equipment.)

#### Flow Cytometry Detection

#### Sample Test Data Acquisition

- 1) Detect samples using this kit (example for Acquire samples

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

带格式的

...

- on a flow cytometer using manual gating analysis).
- 2) Create an FSC-H/SSC-H scatter plot. Adjust the gain of FSC and SSC to gate/identify the cytokine bead population using a rectangular P1 gate, as shown in (Figure 1).




Figure 1

- 3) ~~Establish~~ Create an APC-H/APC-Cy7-H scatter plot. ~~Identify the individual bead populations from the main P1-gated bead cluster population, identify individual bead populations~~ and gate each ~~specific~~-cytokine-specific bead population using rectangular gates, as shown in (Figure 2).




Figure 2

## Data Analysis

- 1) ~~Open the template for detection and s~~Set the data acquisition to collect at least 200 beads per cytokine gate. For example, ~~when analyzing four bead populations, collect if detecting four distinct bead populations,~~ a minimum of  $4 \times 200 = 800$  total beads should be collected. ~~Proceed with data analysis after acquisition is complete.~~
- 2) Calculate the median fluorescence intensity (MFI) of the standard and sample ~~replicate wells,~~ then ~~s~~Subtract the MFI ~~value~~ of the blank well to obtain the corrected MFI ~~value~~.
- 3) ~~Generate a standard curve by plotting concentration on the (x-axis) and MFI on the (y-axis) using a log-log scale, and fit the curve the data with a four~~4-parameter logistic (4-PLK) ~~function model.~~
- 4) If the sample MFI exceeds the upper limit of the standard curve, ~~appropriately~~ dilute the sample ~~appropriately~~ and repeat the measurement. Multiply the calculated concentration by the corresponding dilution factor.
- 5) ~~Data Analysis: Analyze Experimental data can be analyzed using the Elabscience® Aptplex™ Analysis Software or other commercially available compatible data analysis software.~~

4)

## Performance parameters

1. Detection range:

Cytokines	Range of linearity
IL-6	5-5000 pg/mL
IL-1β	5-5000 pg/mL

TNF-α	5-5000 pg/mL
IFN-γ	5-5000 pg/mL
CXCL10 (IP-10)	10-5000 pg/mL
CCL2 (MCP-1)	10-5000 pg/mL
IL-8 (CXCL8)	5-5000 pg/mL
IL-10	5-5000 pg/mL
Cytokines	Range of linearity
IL-17A	5-5000 pg/ml
IL-6	5-5000 pg/ml
TNF-α	5-5000 pg/ml
IFN-γ	5-5000 pg/ml
IL-12p70	5-5000 pg/ml
IL-10	5-5000 pg/ml
CCL5 (RANTES)	10-5000 pg/ml
CXCL10 (IP-10)	10-5000 pg/ml
IL-1β	5-5000 pg/ml
IL-18	10-5000 pg/ml
sCD54 (ICAM-1)	40-10000 pg/ml
IL-2RA (sCD25)	10-5000 pg/ml

2. Limit of Blank (LoB):

带格式的: 上标

带格式的: 上标

2. The LoB for all assays-analyses is  $\leq 8$  pg/mL.
3. Recovery:
3. The mean recovery rate ranges from 70% to 120%.
4. Precision:
4. ~~Both~~ The intra-assay and inter-assay coefficients of variation (CV) are  $\leq 15\%$ .
5. Specificity:
5. No significant cross-reactivity is observed among the different factors-analytes within included in this kit.

## Precautions

- This product is intended for **research use only and must be used** by qualified professionals ~~only~~. Personnel responsible for data interpretation and reporting ~~must receive proper~~ should have appropriate technical training.
- ~~Please observe all safety precautions and adhere to standard laboratory reagent handling procedures~~ Follow standard laboratory safety practices and reagent handling procedures. This product contains fluorescent dyes. Avoid direct contact with skin and eyes, and do not allow ~~prevent~~ contamination of food and beverages, and ~~Always~~ wear appropriate gloves when ~~personal protective equipment, including gloves, during handling~~.
- ~~Incorrect-Improper calibration of the flow cytometer calibration, improper-inadequate compensation for fluorescence compensation-spillover, or imprecise-incorrect gating strategies may lead to erroneous-inaccurate results. Please refer to the instrument's operating manufacturer's manual for and ensure proper calibration procedures to sample acquisition ensure optimal performance prior to use.~~

- ~~Prior to Before~~ use, ~~thoroughly~~ vortex the bead suspension thoroughly to ensure ~~homogeneity-uniform bead dispersion~~ and ~~to~~ prevent bead aggregation, which ~~could lead to uneven bead distribution across wells may affect assay performance~~.
- To prevent cross-contamination, change pipette tips between each well, exercise caution when removing the plate sealer to avoid contact with adjacent wells; ~~always~~ use new/fresh tips ~~when handling for~~ different standards or samples; ~~and avoid generating-bubbles formation~~ during pipetting. ~~The use~~ of a multichannel pipette is recommended for adding wash buffer steps.
- ~~To ensure fluorescence detection quality, protect the reaction from light during all steps involving detection antibodies~~ Protect all reactions involving detection antibodies and SA-PE from light throughout the assay to maintain fluorescence signal integrity.
- Do not mix reagents from different lot numbers. ~~Do not use products from or substitute reagents from other manufacturers interchangeably. Strictly follow the instructions provided in this manual. Use the product only within its stated validity period. Use all components according to the instructions provided in this manual and within their stated expiration dates.~~

带格式的: 缩进: 左侧: 0.74 厘米, 无项目符号或编号

带格式的: 缩进: 左侧: 0.74 厘米, 无项目符号或编号

带格式的: 缩进: 左侧: 0.74 厘米, 无项目符号或编号

带格式的: 缩进: 左侧: 0.74 厘米, 无项目符号或编号