A Reliable Research Partner in Life Science and Medicine

Recombinant Mouse Carbonic Anhydrase II/CA2 Protein (His Tag)

Catalog Number: PKSM040550

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Species Mouse

Source E.coli-derived Mouse Carbonic Anhydrase II/CA2 protein Met 1-Lys 260, with an C-

terminal His

Calculated MW 30.4 kDa Observed MW 30-33 kDa Accession AAH55291.1

Measured by its esterase activity. The specific activity is > 100 pmoles/min/µg. **Bio-activity**

Properties

> 96 % as determined by reducing SDS-PAGE. **Purity** Endotoxin Please contact us for more information.

Storage Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80

°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of

reconstituted samples are stable at < -20°C for 3 months.

This product is provided as lyophilized powder which is shipped with ice packs. Shipping

Lyophilized from sterile 50mM Tris, 150mM NaCl, pH 7.5 **Formulation**


Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants

before lyophilization.

Please refer to the specific buffer information in the printed manual.

Reconstitution Please refer to the printed manual for detailed information.

Data



> 96 % as determined by reducing SDS-PAGE.

Background

Elabscience Bionovation Inc.



A Reliable Research Partner in Life Science and Medicine

The carbonic anhydrases (or carbonate dehydratases) are classified as metalloenzyme for its zinc ion prosthetic group and form a family of enzymes that catalyze the rapid interconversion of carbon dioxide and water to bicarbonate and protons, a reversible reaction that takes part in maintaining acid-base balance in blood and other tissues. The carbonic anhydrasekl (CA) family consists of at least 11 enzymatically active members and a few inactive homologous proteins. Carbonic anhydrase II is one of fourteen forms of human α carbonic anhydrases. Defects in this enzyme are associated with osteopetrosis and renal tubular acidosis. Renal carbonic anhydrase allows the reabsorption of sodium ions in the proximal tubule. Carbonic anhydrase II has been shown to interact with Band 3 and Sodium-hydrogen antiporter 1.

Fax: 1-832-243-6017