The professional cell culture empowers a healthier world

HT-3 Cell Complete Medium

Cat. No. : CM-0630 Size:125mL×4

General Information

Product Form Liquid

Concentration Ready-to-use

Components McCoy's 5A(PM150710)+Nutrients

Bacterial detection Negative

Fungal detection Negative

Mycoplasmal detection Negative

Endotoxin level <3 EU/mL

Shipping Conditions Ice bag

Storage 2-8°C ,Shading Light

Expiry date 6 months

Product Introduction

HT-3 Cell Complete Medium is a ready-to-use medium that can be directly used for the culture of HT-3 cells, which is carefully designed and optimized by Pricella's R&D team based on McCoy's 5A basal medium. This high-quality complete medium has been strictly verified by HT-3 cells which could promote cell proliferation and maintain cell properties.

Guidelines for use

- 1. Pricella's cell culture media undergoes strictquality control to ensure sterility, but may get contaminated during use. Follow these guidelines for sterile handling toavoid contamination.
- 2. Always wipe your gloved hands and work area with 70% ethanol.
- 3. Wipe the outside of the containers, flasks, plates, and dishes with 70% ethanol before placing them in the cell culture hood.
- 4. Use sterile pipette tips and pipettes to work with liquids, and use each pipette tip only once to avoid cross-contamination.Do not unwrap sterile pipettes until they are ready to beused. Keep pipettes and tips within the
- 5. Do not talk while performing sterile procedures and perform your cell culture as efficiently and carefully as possible tominimize contamination.

Quality control

Web: www.pri-cella.com

Standard evaluations for cell culture media are pH, osmolality, endotoxins and sterility testing for liquid products, cell growth experiments .

Matters Needing Attention

- 1. This product is for research use only.
- 2. This product is sterilized by 0.1 μm filtration.
- 3. It is not suitable for long time storage at room temperature.
- 4. This product is a ready-to-use medium. If there is no special need, don't add serum, penicillin and streptomycin. It can be used directly.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Rev V2.7