



The professional cell culture empowers a healthier world

# Ham's F-10 Complete Medium (10)

Cat. No.: PM151110B

Size: 125mL×4

### **General Information**

**Product Form** Liquid Concentration Ready-to-use

Components Ham's F-10[PM151110]+10% Nutrient+1% Supplement

**Bacterial detection** Negative Fungal detection Negative Mycoplasmal detection Negative **Endotoxin level** < 3 EU/mL Shipping Ice bag

2-8°C, Shading Light Storage

**Expiry date** 6 months

## **Product Introduction**

icella Ham's F-10 Nutrient Mixture is designed by Ham in 1963 for serum-free culture of CHO cells. Ham's F-10 is rich in nutrients and is suitable for the culture of human diploid cells, clonal culture of cells under low serum conditions, normal culture of mammalian cells, chromosome analysis of white blood cells, and primary cell culture of rats, rabbits, chickens and other animals.

## **Guidelines for use**

- 1. Pricella's cell culture media undergoes strictquality control to ensure sterility, but may get contaminated during use. Follow these guidelines for sterile handling toavoid contamination.
- 2. Always wipe your gloved hands and work area with 70% ethanol.
- 3. Wipe the outside of the containers, flasks, plates, and dishes with 70% ethanol before placing them in the cell culture hood.
- 4. Use sterile pipette tips and pipettes to work with liquids, and use each pipette tip only once to avoid cross contamination. Do not unwrap sterile pipettes until they are ready to beused. Keep pipettes and tips within the clean work area.
- 5. Do not talk while performing sterile procedures and perform your cell culture as efficiently and carefully as possible tominimize contamination.

# **Quality control**

Standard evaluations for cell culture media are pH, osmolality, endotoxins and sterility testing for liquid products, cell growth experiments.

### **Notes**

- 1. This product is for research use only.
- 2. This product is sterilized by 0.1 µm filtration.
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.
- 4. It is not suitable for long time storage at room temperature.
- 5. This product is a ready-to-use medium. If there is no special need, don't add serum, penicillin and streptomycin. It can be used directly.

Tel:1-832-243-6086 Toll-free:1-888-852-8623 Fax:1-832-243-6017 Email: techsupport@pri-cella.com Web: www.pri-cella.com Rev. V2.9