The professional cell culture empowers a healthier world

Low Calcium Delklund's Solution

Cat. No.: PB180534

Size: 500mL

General Information

Product Form Liquid Concentration 1 × pН 7.2 - 7.4NaCl 9000 mg/L **KCI** 420 mg/L CaCl₂ 30 mg/L Sodium Bicarbonate (NaHCO₃) 500 mg/L **D-Glucose** 500 mg/L Storage 2-30°C

rocell system **Shipping** Room Temperature

Expiration date 24 months

Background

Dejalon's solution, also known as rat uterine fluid, is commonly used for cleaning and infiltrating rat uterine tissue.

Guidelines for use

- Buffered saline solutions have many uses in cell culture protocols such as cell washing, diluents, or as research sample holding solution.
- 2. Live cellular suspensions prepared in a buffered saline solution should not be stored longer than several hours, as cellular viability may decrease. Therefore, the suitability of the buffered saline solution for a specific cell type should be tested before use.

Notes

- This product is for research use only.
- This product is sterilized by 0.1 µm filtration.
- It is necessary to pay attention to the aseptic operation and avoid the contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Web: www.procellsystem.com Email: techsupport@procellsystem.com