The professional cell culture empowers a healthier world

## **Smooth Muscle Cell Medium (SMCM)**

Cat. No.: CM-ZY006

Size: 500mL

#### **General Information**

Product Form Liquid
Concentration 1×
Size 500mL

The proportion of complete 92% Basic Culture Medium for Smooth Muscle Cells + 5% Nutrient for medium preparation

Smooth Muscle Cells (for Human and Rat) + 2% Supplement for Smooth

Muscle Cells (for Human and Rat) + 1% Additive for Smooth Muscle Cells

**Fungal detection** Negative Negative

Cells growth experiment The cells grow well and have a normal morphology

Negative

Endotoxin level < 3 EU/mL

Storage Store at the temperature indicated on the label, shading light

Shipping Ice bag

## **Background**

Mycoplasmal detection

SMCM has been meticulously optimized by the Procell's technical team. After long-term testing, this product can maintain the best growth state of Smooth Muscle Cells (for Human and Rat). This product already contains all the components required for the growth of Smooth Muscle Cells (for Human and Rat). If added in the recommended proportion, no other components need to be added and it can be directly used for the in vitro culture of Smooth Muscle Cells (for Human and Rat). This product is for further research use only and must not be used for diagnosis, treatment, clinical, household or any other purpose.

# Components of the culture medium

Components	Volume	Storage	Expiration date
Basic Culture Medium for Smooth Muscle Cells	500mL	2-8°C, shading light	12 months
Nutrient for Smooth Muscle Cells (for Human and Rat)	25mL	-5~-20°C, shading	60 months
Supplement for Smooth Muscle Cells (for Human and	10mL	-5~-20°C, shading	12 months
Additive for Smooth Muscle Cells	5mL	-5~-20°C, shading	12 months

Supplement for Smooth Muscle Cells (for Human and Rat), containing(part):	Volume	Final Concentration in Complete Medium	Storage
EGF	50μL	1 ng/mL	-5~-20°C, shading light
FGF basic	100µL	2 ng/mL	-5~-20°C, shading light
Insulin	250µL	5 μg/mL	-5~-20°C, shading light

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Web: www.procellsystem.com Email: techsupport@procellsystem.com





The professional cell culture empowers a healthier world

### Note

- 1. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.
- 2. Nutrient for Smooth Muscle Cells (for Human and Rat), Supplement for Smooth Muscle Cells (for Human and Rat) and Additive for Smooth Muscle Cells should be stored at the temperature indicated on the label. It can be used after dissolving at 2-8°C. To avoid repeated freezing and thawing, it is recommended to store it in small quantities after subpackaging.
- 3. This product is packaged in a single component form. It needs to be prepared into a complete culture before use. The validity period of the complete culture medium is 3 months.
- 4. In order to maintain the best use effect of this product, it should not be placed in room temperature or high temperature environment for a long time.
- 5. All products should be used within their expiration dates. Once the expiration date is exceeded, they must be discarded.
- 6. This product is only for scientific research or further research use, not for diagnosis and treatment.







Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Web: <u>www.procellsystem.com</u> Email: <u>techsupport@procellsystem.com</u>