



A Reliable Research Partner in Life Science and Medicine

Elab Fluor[®] Violet 450 Rat IgM, κ Isotype Control[RTK2118]

Catalog Number: E-AB-F09772Q

Note: Centrifuge before opening to ensure complete recovery of vial contents.

| _ | | 4.0 | |
|----|-------|------|-----|
| | escri | MTI. | on. |
| ш, | COUL | | UII |

Host Rat Isotype Rat IgM, ĸ RTK2118 Clone No.

Conjugation Elab Fluor® Violet 450

Conjugation Information Elab Fluor® Violet 450 is designed to be excited by the violet laser (405 nm) and

detected using an optical filter centered near 450 nm (e.g., a 450/45 nm bandpass filter).

Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein Storage Buffer

protectant.

Applications Recommended usage

FCM Each lot of this antibody is quality control tested by flow cytometric analysis. The amount

> of the reagent is suggested to be used 5 μ L of antibody per test (million cells in 100 μL staining volume or per 100 μL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for

individual use.

Preparation & Storage

Keep as concentrated solution. Storage

This product can be stored at 2-8°C for 12 months. Please protected from prolonged

exposure to light and do not freeze.

Shipping Ice bag

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017 Web:www.elabscience.com