Recombinant Mouse ICAM-2/CD102 Protein (Fc Tag)

Catalog Number: PKSM041208



Note: Centrifuge before opening to ensure complete recovery of vial contents.

| Description | |
|-------------|--------------|
| Species | Mouse |
| Mol_Mass | 49.7 kDa |
| Accession | P35330 |
| | 3.7 . 11.1 . |

Bio-activity Not validated for activity

| т. | | | | | , . | |
|----|----|---|----|----------------|-----|----|
| 12 | rn | n | • | T ^a | т | es |
| | | | w. | | ш | |

Description

> 95 % as determined by reducing SDS-PAGE. **Purity**

< 1.0 EU per ug of the protein as determined by the LAL method. Endotoxin

Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 Storage

°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of

reconstituted samples are stable at < -20°C for 3 months.


This product is provided as lyophilized powder which is shipped with ice packs. Shipping Lyophilized from a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0. Formulation Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants

before lyophilization.

Please refer to the specific buffer information in the printed manual.

Please refer to the printed manual for detailed information. Reconstitution

Data



Background

ICAM-2 is a 55-65 kD transmembrane glycoprotein possessing 2 extracellular Ig domains, a single transmembrane domain, and a short 26-amino acid cytoplasmic domain. ICAM-2 is expressed on most leukocytes, and is strongly expressed on vascular endothelial cells. Interactions of ICAM-2 and the integrin receptors mediate cell adhesion in a wide range of lymphocyte, monocyte, natural killer cell, and granulocytewith other cells, and play important roles in many adhesion-dependent immune and inflammation responses, such as T cell aggregation, NK-cell cytotoxicity and migration, lymphocyte recirculation, etc. Serum levels of ICAM-2 correlated significantly with the inflammatory and course sequences of trichinosis in mice and had a similar relation with blood eosinophilia. So, estimation of ICAM-2 serum levels may prove useful in diagnosis of trichinosis recent infections, and in monitoring the prognosis and response to treatment.

For Research Use Only