

### **HT-1197 Cell Complete Medium**

Cat. No.: CM-0824 Size:125mL×4

#### **General Information**

**Product Form** Liquid

Concentration Ready-to-use

Components MEM, with NEAA(PM150410)+Nutrients+1 mM Sodium Pyruvate (PB180422)

**Bacterial detection** Negative

**Fungal detection** Negative

Mycoplasmal detection Negative

bsci€3 EÛ/mL

**Shipping Conditions** 

**Endotoxin level** 

Ice bag

Storage

2-8°C, Shading Light

**Expiry date** 6 months

#### **Product Introduction**

HT-1197 Cell Complete Medium is a ready-to-use medium that can be directly used for the culture of HT-1197 cells, which is carefully designed and optimized by Pricella's R&D team based on MEM, with NEAA basal medium. This high-quality complete medium has been strictly verified by HT-1197 cells which could promote cell proliferation and maintain cell properties.

### **Guidelines for use**

- Pricella's cell culture media undergoes strictquality control to ensure sterility, but may get contaminated during use. Follow these guidelines for sterile handling to avoid contamination.
- Always wipe your gloved hands and work area with 70% ethanol.
- Wipe the outside of the containers, flasks, plates, and dishes with 70% ethanol before placing them in the cell culture hood.
- Use sterile pipette tips and pipettes to work with liquids, and use each pipette tip only once to avoid cross-contamination. Do not unwrap sterile pipettes until they are ready to beused. Keep pipettes and tips within the clean work area.
- Do not talk while performing sterile procedures and perform your cell culture as efficiently and carefully as possible tominimize contamination.

# **Quality control**

Standard evaluations for cell culture media are pH, osmolality, endotoxins and sterility testing for liquid products, cell growth experiments.

# **Matters Needing Attention**

- 1. This product is for research use only.
- 2. This product is sterilized by 0.1 µm filtration.
- 3. It is not suitable for long time storage at room temperature.
- This product is a ready-to-use medium. If there is no special need, don't add serum, penicillin and streptomycin. It can be used directly.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017 Web: www.pri-cella.com

Rev V2.8

Email: techsupport@pri-cella.com