

Leibovitz's L-15, with HEPES

Cat. No.: PM151015

Size: 500mL

General Information

Product Form Liquid Concentration 1 × 7.2 - 7.4pН Negative Sodium Bicarbonate (NaHCO3) Negative L-Glutamine 300 mg/L**HEPES** 5958 mg/L Phenol red 10 mg/L**D**-galactose 900 mg/L Sodium pyruvate 550 mg/L **Balanced salt** Phosphate

ricella Elabscience Storage 2-8°C, Shading Light **Shipping** Room Temperature

Expiration date 24 months

Background

MEM medium (Minimum Essential Medium) was developed on the basis of Eagle basic medium. It is one of the most basic and widely used culture medium, and one of the most commonly used culture medium in animal cell culture. MEM medium contains 12 kinds of essential amino acids, glutamine and 8 vitamins, which is simple, mainly used in the culture of adherent cells.

MEM (contain NEAA) medium is added L-alanine, L-glutamic acid, L-asparagine, L-aspartic acid, L-proline, L-serine and glycine on the basis of MEM medium. These 7 kinds of NEAA, can reduce the side effects of producing non-essential amino acids during cell culture and promote cell proliferation and metabolism.

This product contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors. Therefore, the product should be used by Elabscien with serum or serum-free additives.

Guidelines for use

Pricella's cell culture media undergoes strict quality control to ensure sterility, but may get contaminated during use. Follow these guidelines for sterile handling to avoid contamination.

- Always wipe your gloved hands and work area with 70% ethanol.
- Wipe the outside of the containers, flasks, plates, and dishes with 70% ethanol before placing them in the cell culture hood.
- Use sterile pipette tips and pipettes to work with liquids, and use each pipette tip only once to avoid

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017 Web: www.pri-cella.com Email: techsupport@pri-cella.com





cross-contamination. Do not unwrap sterile pipettes until they are ready to be used. Keep pipettes and tips within the clean work area.

4. Do not talk while performing sterile procedures and perform your cell culture as efficiently and carefully as possible to minimize contamination.

Quality control

Standard evaluations for cell culture media are pH, osmolality, endotoxins and sterility testing for liquid products, cell growth experiments .

Notes

- 1. This product is for research use only.
- 2. This product is sterilized by 0.1 µm filtration.
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.









Rev V1.1