

PVDF Membrane (0.22 μ m)

Cat. No: E-BC-R329

Size: 10 pieces / 20 pieces



Product Content

Cat	Products	Size		Storage
E-BC-R329	PVDF Membrane (0.22 μ m, 8.5 cm $\times$ 6 cm)	10 pieces	20 pieces	RT
	PVDF Membrane (0.22 μ m, 13.25 cm $\times$ 15 cm)	10 pieces	20 pieces	

Binding ability: Insulin: 262 μ g/cm²

BSA: 340 μ g/cm²

Goat IgG: 448 μ g/cm²

Introduction

Compared with nitrocellulose membranes, hydrophobic polyvinylidene fluoride (PVDF) membranes have many advantages, for example, good mechanical strength, strong protein retention, low background, extensive solvent compatibility, and excellent coloring ability. The PVDF membrane with a pore size of 0.22 μ m can be used for most western blots and is recommended for proteins with a molecular weight less than 20 kDa.

Storage

Store at room temperature for 12 months.

Cautions

1. This kit is for research use only.
2. Please wear gloves to prevent contamination of the PVDF membrane.
3. During the western blot experiment, activate the PVDF membrane by soaking in methanol for 1 min in advance.
4. Do not dry the PVDF membrane during the experiment.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Rev. V1.3