

## Vero Serum-Free Medium

**Cat. No. : SF4001**

**Size: 500mL**

### Description

Vero Serum-Free Medium is a chemically defined medium designed and developed specifically for Vero cells. It is suitable for Vero cell adherent culture. Vero Serum-Free Medium contains L-Glutamine, which generally does not require additional addition.

### General Information

Product Form	Liquid
Sterile Filtered	Yes
Specification	500mL®
pH	7.2-7.6
Storage	2-8°C, Shading Light
Shipping	Ice bag
Expiration date	12 months

### Use Instructions

This product can be used directly for serum-free cultured Vero cells. If the culture is transferred from serum-containing medium to a serum-free medium, it is necessary to domesticate Vero cells using a serum-free medium. The taming gradient is 2%-1%-0.5%-0% (the proportion of serum added to Vero Serum-Free Medium), and serum-free cell domestication can be completed within 1-2 weeks.

### Notes

1. This product has been filtered and sterilized. Pay attention to aseptic operation to avoid contamination.
2. To maintain the best use of this product, do not perform freeze-thaw treatment.
3. This product is for research use only, and cannot be used for clinical diagnosis or treatment.