

The professional cell culture empowers a healthier world

TNM-FH medium for insect Cells

Cat. No.: PM152010

Size: 500mL

General Information

Product Form Liquid Concentration 1 × 6.1-6.4 pН **D-Glucose** 700 mg/L₂ **HEPES** Negative 600 mg/L L-Glutamine NaHCO₃ 350 mg/L Phenol red Negative

Storage 2-8°C, Shading Light procell system Shipping **Room Temperature**

Expiration date 12 months

Background

TNM-FH medium is an improved Grace's insect medium. After proper supplementation, it can be used for the culture of various Lepidoptera insect cells.

TNM-FH medium contains amino acids, vitamins, inorganic salts and other ingredients required for various types of cell culture, but does not contain protein, lipids or any growth factors, so this product needs to be used with serum or serum-free additives.

Guidelines for use

Procell's cell culture media undergoes strict quality control to ensure sterility, but may get contaminated during use. Follow these guidelines for sterile handling to avoid contamination.

- Always wipe your gloved hands and work area with 70% ethanol.
- 2. Wipe the outside of the containers, flasks, plates, and dishes with 70% ethanol before placing them in the cell culture hood.
- 3. Use sterile pipette tips and pipettes to work with liquids, and use each pipette tip only once to avoid cross-contamination. Do not unwrap sterile pipettes until they are ready to be used. Keep pipettes and tips within the clean work area.
- 4. Do not talk while performing sterile procedures and perform your cell culture as efficiently and carefully as possible to minimize contamination.

Quality control

Standard evaluations for cell culture media are pH, osmolality, endotoxins and sterility testing for liquid products, cell growth experiments.

Notes

- This product is for research use only.
- 2. This product is sterilized by 0.1 µm filtration.
- It is necessary to pay attention to the aseptic operation and avoid the contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Web: www.procellsystem.com Email: techsupport@procellsystem.com