

Smooth Muscle Cell Medium (SMCM)

Cat. No.: CM-ZY006 Size: 500mL

General Information

Product Form	Liquid
Concentration	1×
Size	500mL
The proportion of complete	e 92% Basic Culture Medium for Smooth Muscle Cells + 5% Nutrient for Smooth
medium preparation	Muscle Cells (for Human and Rat) + 2% Supplement for Smooth Muscle Cells
ice	(for Human and Rat) + 1% Additive for Smooth Muscle Cells
Bacterial detection	Negative
Fungal detection abscience	Negative
Mycoplasmal detection	Negative
Cells growth experiment	The cells grow well and have a normal morphology
Endotoxin level	< 3 EU/mL 11 0 ®
Storage	Store at the temperature indicated on the label, shading light
Shipping	Ice bag
Background	by Elabse

Background

SMCM has been meticulously optimized by the Pricella's technical team. After long-term testing, this product can maintain the best growth state of Smooth Muscle Cells (for Human and Rat). This product already contains all the components required for the growth of Smooth Muscle Cells (for Human and Rat). If added in the recommended proportion, no other components need to be added and it can be directly used for the in vitro culture of Smooth Muscle Cells (for Human and Rat). This product is for further research use only and must not be used for diagnosis, treatment, clinical, household or any other purpose.

Components of the culture medium

Components		Volume	Storage	Expiration date
Basic Culture Medium for Smooth Muscle Cells			2-8°C, shading light	12 months
Nutrient for Smooth Muscle Cells (for Human and Rat)			-5~-20°C, shading	60 months
Supplement for Smooth Muscle Cells (for Human and Rat)			-5~-20°C, shading	12 months
Additive for Smooth muscle cells			-5~-20°C, shading	12 months
		PLI	science	
Supplement for Smooth Muscle Cells (for Human and Rat), containing(part):	Volume	Final Conce Complete M	ntration in Stor	age
EGF	50µL	1 ng/mL	-5~-2	20°C, shading light
FGF basic	100µL	2 ng/mL	-5~-2	20°C, shading light
Insulin	250µL	5 μg/mL	-5~-2	20°C, shading light



- 1. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.
- 2. Nutrient for Smooth Muscle Cells (for Human and Rat), Supplement for Smooth Muscle Cells (for Human and Rat) and Additive for Smooth Muscle Cells should be stored at the temperature indicated on the label. It can be used after dissolving at 2-8°C. To avoid repeated freezing and thawing, it is recommended to store it in small quantities after subpackaging.
- 3. This product is packaged in a single component form. It needs to be prepared into a complete culture before use. The validity period of the complete culture medium is 3 months.
- 4. In order to maintain the best use effect of this product, it should not be placed in room temperature or high temperature environment for a long time.
- 5. All products should be used within their expiration dates. Once the expiration date is exceeded, they must be discarded.
- 6. This product is only for scientific research or further research use, not for diagnosis and treatment.

Pricella[®] by Elabscience



