

MEM, with NEAA, P/S, without L-glutamine

Cat. No: PM150412A Size: 500mL

General Information

Product Form	Liquid
Concentration	1×
рН	7.2-7.4
D-Glucose	1000mg/L
HEPES Elabscience	Negative
L-Glutamine	Negative
NaHCO3	2200mg/L
Phenol red	10mg/L
Penicillin G	100kU/L
Earle's Balanced Salt	Positive
NEAA	Positive
Streptomycin Sulfate	100mg/L
Storage	2-8 $^\circ\!\!\mathbb{C}$,Shading Light
Shipping	RT
Expiration date	24 months



Background

MEM medium (Minimum Essential Medium) was developed on the basis of Eagle basic medium. It is one of the most basic and widely used culture medium, and one of the most commonly used culture medium in animal cell culture. MEM medium contains 12 kinds of essential amino acids, glutamine and 8 vitamins, which is simple, mainly used in the culture of adherent cells.

MEM (contain NEAA) medium is added L-alanine, L-glutamic acid, L-asparagine, L-aspartic acid, L-proline, Lserine and glycine on the basis of MEM medium. These 7 kinds of NEAA, can reduce the side effects of producing non-essential amino acids during cell culture and promote cell proliferation and metabolism. This product contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors. Therefore, the product should be used by Elabscience with serum or serum-free additives.

Notes

- 1. This product is for research use only.;
- 2. This product is sterilized by 0.1µm filtration.;
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.