

he professional cell culture enpowers a healthier world

## **Advanced Tyrode's Solution, with HEPES**

Cat. No: PB180340

Size: 500mL

## **General Information**

Product Form

Concentration

pH

7.2-7.4

D-Glucose

HEPES

5 mM

NaHCO3

Negative

Phenol red

Negative

CaCl2 (anhydrous)

Liquid

1×

7.2-7.4

Negative

200mg/L

Storage 2-30°C ,Shading Light

214mg/L

Shipping RT

**Expiration date** 24 months



## **Background**

MgCl2•6H2O

Balanced Salt Solution (Physiological Solution) have the properties of buffer solution (regulate pH), normal saline (maintain osmotic pressure) and culture medium (provide nutrition). It can meet the basic needs of survival and metabolism of tissues, organs or cells in vitro. A small amount of phenolic red was added to some equilibrium salt solutions to indicate the pH change of the solution.

Tyrode's solution is a kind of balanced salt solution in biochemistry. The main components are NaCl, KCl, NaH2PO4, MgCl2, CaCl2, NaHCO3 and glucose. Tyrode's solution is often used for perfusion and tissue cleaning in the body of mammals, and is especially suitable for isolated constant-muscle experiments to maintain the normal physiological function of isolated enteromuscles.

Hepes was added to advanced Tyrode's solution and the buffer capacity of the solution was greatly improved.

## Notesby

- 1. This product is for research use only.;
- 2. This product is sterilized by 0.1µm filtration.;
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Web: www.pri-cella.com Email: techsupport@pri-cella.com