

The professional cell culture empowers a healthier world

rocellsystem

## **Ringer's Solution, for amphibians**

Cat. No.: PB180342

Size: 500mL

## **General Information**

**Product Form** Liquid Concentration 1 × 7.2-7.4 pН 2000 mg/L **D-Glucose** Negative **HEPES** 200 mg/L NaHCO<sub>3</sub> Phenol red Negative CaCl<sub>2</sub> (anhydrous) 120 mg/L 2-30°C **Storage** 

Shipping Room Temperature

**Expiration date** 36 months

## **Background**

Balanced Salt Solution (Physiological Solution) have the properties of buffer solution (regulate pH), normal saline (maintain osmotic pressure) and culture medium (provide nutrition). It can meet the basic needs of survival and metabolism of tissues, organs or cells in vitro. A small amount of phenolic red was added to some equilibrium salt solutions to indicate the pH change of the solution. Ringer's Solution was invented by British physiologist Ringer. It is a kind of balanced salt solution, which is close to the amphibian's internal environment. In vitro experiments, it can be used to prolong the beating time of frog's heart in vitro and maintain the physiological activity of other isolated tissues and organs in vitro.

## **Notes**

This product is for research use only.

Web: www.procellsystem.com

- 2. This product is sterilized by 0.1 µm filtration.
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Email: <a href="mailto:techsupport@procellsystem.com">techsupport@procellsystem.com</a>