

The professional cell culture empowers a healthier world

Procell system

ell systen

## **Tyrode's Solution**

Cat. No.: PB180338

Size: 500mL

## **General Information**

**Product Form** Liquid Concentration 1 × pН 7.2 - 7.4**D-Glucose** 1000 mg/L **HEPES** Negative 1000 mg/L NaHCO<sub>3</sub> Phenol red Negative

CaCl<sub>2</sub> (anhydrous) 200 mg/L MgCl<sub>2</sub>•6H<sub>2</sub>O 214 mg/L 2-30°C **Storage** 

Room Temperature **Shipping** 

**Expiration date** 36 months

## **Background**

Balanced Salt Solution (Physiological Solution) have the properties of buffer solution (regulate pH), normal saline (maintain osmotic pressure) and culture medium (provide nutrition). It can meet the basic needs of survival and metabolism of tissues, organs or cells in vitro. A small amount of phenolic red was added to some equilibrium salt solutions to indicate the pH change of the solution. Tyrode's solution is a kind of balanced salt solution in biochemistry.

The main components are NaCl, KCl, NaH<sub>2</sub>PO<sub>4</sub>, MqCl<sub>2</sub>, CaCl<sub>2</sub>, NaHCO<sub>3</sub> and glucose. Tyrode's solution is often used for perfusion and tissue cleaning in the body of mammals, and is especially suitable for isolated constant-muscle experiments to maintain the normal physiological function of isolated enteromuscles.

## **Notes**

- This product is for research use only.
- This product is sterilized by 0.1 µm filtration.
- It is necessary to pay attention to the aseptic operation and avoid the contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017 Web: www.procellsystem.com

Email: techsupport@procellsystem.com