

## Smooth Muscle Cell Medium-Mouse (SMCM-m)

Cat. No. : CM-ZY007

Size: 500mL

### General Information

<b>Product Form</b>	Liquid
<b>Concentration</b>	1 ×
<b>Size</b>	500mL
<b>The proportion of complete medium preparation</b>	87% Basic Culture Medium for Smooth Muscle Cells + 10% Nutrient for Smooth Muscle Cells (for Mouse) + 2% Supplement for Smooth Muscle Cells (for Mouse) + 1% Additive for Smooth Muscle Cells
<b>Bacterial detection</b>	Negative
<b>Fungal detection</b>	Negative
<b>Mycoplasma detection</b>	Negative
<b>Cells growth experiment</b>	The cells grow well and have a normal morphology
<b>Endotoxin level</b>	< 3 EU/mL
<b>Storage</b>	Store at the temperature indicated on the label, shading light
<b>Shipping</b>	Ice bag

### Background

SMCM-m has been meticulously optimized by the Pricella's technical team. After long-term testing, this product can maintain the best growth state of Smooth Muscle Cells (for Mouse). This product already contains all the components required for the growth of Smooth Muscle Cells (for Mouse). If added in the recommended proportion, no other components need to be added and it can be directly used for the in vitro culture of Smooth Muscle Cells (for Mouse). This product is for further research use only and must not be used for diagnosis, treatment, clinical, household or any other purpose.

### Components of the culture medium

Components	Volume	Storage	Expiration date
Basic Culture Medium for Smooth Muscle Cells	500mL	2-8°C, shading light	12 months
Nutrient for Smooth Muscle Cells (for Mouse)	50mL	-5~-20°C, shading light	60 months
Supplement for Smooth Muscle Cells (for Mouse)	10mL	-5~-20°C, shading light	12 months
Additive for Smooth Muscle Cells	5mL	-5~-20°C, shading light	12 months

Supplement for Smooth Muscle Cells (for Mouse), containing(part):	Volume	Final Concentration in Complete Medium	Storage
FGF basic	100μL	2 ng/mL	-5~-20°C, shading light
Insulin	250μL	5 μg/mL	-5~-20°C, shading light

### Note

1. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.
2. Nutrient for Mouse Smooth Muscle Cells, Supplement for Mouse Smooth Muscle Cells and Additive for

Smooth Muscle Cells should be stored at the temperature indicated on the label. It can be used after dissolving at 2-8°C. To avoid repeated freezing and thawing, it is recommended to store it in small quantities after subpackaging.

3. This product is packaged in a single component form. It needs to be prepared into a complete culture before use. The validity period of the complete culture medium is 3 months.
4. In order to maintain the best use effect of this product, it should not be placed in room temperature or high temperature environment for a long time.
5. All products should be used within their expiration dates. Once the expiration date is exceeded, they must be discarded.
6. This product is only for scientific research or further research use, not for diagnosis and treatment.

Pricella®  
by Elabscience

Pricella®  
by Elabscience

Pricella®  
by Elabscience

Pricella®  
by Elabscience