

Grace's supplement medium for insect Cells

Cat. No.: PM152012

Size: 500mL

General Information

Product Form Liquid Concentration $1 \times$ 6.1-6.4 рH 700 mg/L **D-Glucose Sodium Bicarbonate (NaHCO3)** 350 mg/L L-Glutamine Flabs 600 mg/L **HEPES** Negative D(-)-Fructose 400 mg/L ricella D(+)-Sucrose 26680 mg/L **Hydrolysed Protein** 3303 mg/L 2-8°C, Shading Light Flabscience Yeast

Storage

Shipping Room Temperature

Expiration date 12 months

Background

Grace's supplement medium for insect cells was originally designed for cell growth of Antherea eucalypti, and is widely used for the cell growth of Spodoptera frugiperda, such as Sf9 and Sf21 cells. Sf9 and Sf21 cells are commonly used in baculovirus expression vector system (BEVS) to express recombinant proteins.

This product contains hydrolyzed milk protein, yeast powder, L-glutamine and sodium bicarbonate. In addition, this product needs to be used with serum or serum-free additives.

Pricella's cell culture media undergoes strict quality control to ensure sterility, but may get contaminated during use. Follow these guidelines for sterile handling to avoid contamination.

- Always wipe your gloved hands and work area with 70% ethanol.
- Wipe the outside of the containers, flasks, plates, and dishes with 70% ethanol before placing them in the 1. cell culture hood.
- Use sterile pipette tips and pipettes to work with liquids, and use each pipette tip only once to avoid crosscontamination. Do not unwrap sterile pipettes until they are ready to be used. Keep pipettes and tips within
- Do not talk while performing sterile procedures and perform your cell culture as efficiently and carefully as possible to minimize contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Email: techsupport@pri-cella.com Web: www.pri-cella.com



Quality control

Standard evaluations for cell culture media are pH, osmolality, endotoxins and sterility testing for liquid products, cell growth experiments .

Notes

- 1. This product is for research use only.
- 2. This product is sterilized by 0.1 µm filtration.
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.









 Toll-free: 1-888-852-8623
 Tel: 1-832-243-6086
 Fax: 1-832-243-6017

 Web: www.pri-cella.com
 Email: techsupport@pri-cella.com