

786-O (3D Spheroid-Specific) Cell Complete Medium

Cat. No.: CM-0010-3D

Size: 125mL×4

General Information

Product Form	Liquid
Concentration	Ready-to-use
Components	Complete culture medium for 786-O (3D)
Bacterial detection	Negative
Fungal detection	Negative
Mycoplasmal detection	Negative
Endotoxin level	< 3 EU/mL
Shipping	Ice bag
Storage	-5~-20°C, shading light, After thawing, store at 2-8°C for up to 1 month
Expiration date	6 months

Product Introduction

The 786-O (3D Spheroid-Specific) Cell Complete Medium is an optimized, serum-free, sterile, ready-to-use medium system. It is specifically developed to support the efficient formation, growth, and maintenance of three-dimensional (3D) tumor spheroids from 786-O human renal cancer cells under suspension conditions.

By precisely mimicking the microenvironment of cancer stem cells (CSCs), this product effectively enriches and selects 786-O CSC populations that possess self-renewal capacity and high tumorigenic potential. It serves as an ideal tool for studying renal cancer stem cell biology, drug resistance, tumor recurrence, and metastatic mechanisms, providing a reliable in vivo-like model for your preclinical research.

Notes

1. This product is only used for scientific research or further research, not for diagnosis and treatment.
2. All operations should be performed in a biological safety cabinet. Wear sterile gloves, a mask, and a lab coat. Maintain aseptic technique to avoid contamination.
3. Long-term storage may result in slight precipitation. Gently invert the bottle to redissolve before use.
4. Use the product only within the indicated shelf life. Discard any expired product.
5. Avoid repeated freeze-thaw cycles.