

Neurobasal-A medium, powder

Cat. No.: PM151236P Size: $5\times1L/1\times10L/1\times50L$

General Information

Product Form Powder **D-Glucose** 4500 mg/L 12.92 g,/L Concentration 2600 mg/L by Elabscienc negative L-Glutamine NaHCO₃ Phenol red 8.1 mg/L Sodium pyruvate negative

Storage 2-8°C, Shading Light **Shipping** Room Temperature

Expiration date 36 months



Neurobasal-A medium is a basic neuron culture medium that does not contain L-glutamine and requires the addition of B-27 neuron Supplement with vitamin A, 50× (product number: PB180637). It is suitable for short-term or long-term maintenance and cultivation of neuronal cells in newborn and mature mice.

Pricella by Elabscience

Preparation method

- The preparation water should be purified water, ultra-pure water or water for injection (WFI), and the water temperature should be controlled between 20-30°C during the preparation process.
- Measure 90% of the final volume preparation water to the solution preparation system. Start stirring, and avoid generating bubbles. For example, if 1 L is required, add 900 mL of preparation water here. And it's recommended that the power output per unit volume (P/V) of the mixing system is greater than 10 W/m^3 .
- Weigh the appropriate amount of powder according to the concentration of 12.92 g/L accurately, and add it to the container prepared in step 2. Stir for more than 20 minutes dissolve all powder completely.
- After the solution is clear, add NaHCO3 at a concentration of 2.2 g/L, continue stirring for 5-10 minutes until dissolved, then add ultra pure water to adjust the volume to the 100% of required.
- If necessary, adjust the pH to 7.10-7.40 with 1 mol/L NaOH solution or 1 mol/L HCl solution. Since filtration will slightly increase the pH, the pH value here is lower than the target pH value (7.10-7.50).
- The prepared solution should be sterilized using a 0.2 µm pore size filter membrane under positive pressure (ensure aseptic technique).
- After filtration, a small amount of liquid culture medium can betaken for quality inspection, and use only after passing the test.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017 Web: www.pri-cella.com Email: techsupport@pri-cella.com

Rev V2.4





The filtered liquid medium should be used immediately or stored in glass bottles, culture medium bottles (PET), or single-use storage bags with an oxygen-barrier coating at 2-8°C away from light. The liquid medium has a shelf life of 1 year under these conditions.

Note

- 1. Please wear a lab coat and use disposable gloves and a mask during operation.
- 2. To ensure the optimal performance of this product, please strictly adhere to the recommended storage conditions for its preservation.
- 3. This product is intended for scientific research exclusively or as a raw material in the production process, and must not be applied for clinical diagnosis or treatment.









Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017 Web: www.pri-cella.com

Email: techsupport@pri-cella.com