

### Trypsin Solution Ingredient Formula

| Ingredient (mg/L) / Cat. No.                         | PB180219 | PB180220 | PB180221 | PB180222 | PB180223 | PB180224 |
|------------------------------------------------------|----------|----------|----------|----------|----------|----------|
| <b>Inorganic Salt</b>                                |          |          |          |          |          |          |
| KCl                                                  | 201      | 400      | 400      | 400      | 400      | 400      |
| KH <sub>2</sub> PO <sub>4</sub>                      | 272      | 60       | 60       | 60       | 60       | 60       |
| NaHCO <sub>3</sub>                                   | —        | 350      | 350      | 350      | 350      | 350      |
| NaCl                                                 | 8006     | 8000     | 8000     | 8000     | 8000     | 8000     |
| Na <sub>2</sub> HPO <sub>4</sub> •12H <sub>2</sub> O | 2865     | 121      | 121      | 121      | 121      | 121      |
| <b>Others</b>                                        |          |          |          |          |          |          |
| Trypsin (1:250)                                      | 500      | 500      | 500      | 500      | 2500     | 2500     |
| EDTA•2Na•2H <sub>2</sub> O                           | 200      | —        | 200      | 200      | —        | —        |
| D-Glucose                                            | —        | 1000     | 1000     | 1000     | 1000     | 1000     |
| Phenol Red                                           | —        | —        | —        | 10       | —        | 10       |

| Ingredient (mg/L) / Cat. No.                         | PB180225 | PB180226 | PB180227 | PB180228 | PB180320 |
|------------------------------------------------------|----------|----------|----------|----------|----------|
| KCl                                                  | 201      | 400      | 400      | 201      | 201      |
| KH <sub>2</sub> PO <sub>4</sub>                      | 272      | 60       | 60       | 272      | 272      |
| NaHCO <sub>3</sub>                                   | —        | 350      | 350      | —        | —        |
| NaCl                                                 | 8006     | 8000     | 8000     | 8006     | 8006     |
| Na <sub>2</sub> HPO <sub>4</sub> •12H <sub>2</sub> O | 2865     | 121      | 121      | 2865     | 2865     |
| Trypsin (1:250)                                      | 2500     | 2500     | 2500     | 2500     | —        |
| EDTA•2Na•2H <sub>2</sub> O                           | 200      | 200      | —        | —        | 200      |
| D-Glucose                                            | —        | 1000     | 1000     | —        | —        |
| Phenol Red                                           | —        | 10       | —        | —        | —        |