



A Reliable Research Partner in Life Science and Medicine

Human ApoD Antibody Pair Set

Catalog No. E-KAB-0506 Applications ELISA

Synonyms APOD

Kit components & Storage

Title	Specifications	Storage
Human ApoD Capture Antibody	1 vial, 100 μ g	Store at -20℃ for one year. Avoid
		freeze/thaw cycles.
Human ApoD Detection Antibody	1 vial, 50 μL	Store at -20°C for one year. Avoid
(Biotin)		freeze/thaw cycles.

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Items		Characteristic (E-KAB-0506)		
		Human ApoD Capture Antibody	Human ApoD Detection Antibody	
			(Biotin)	
Immunogen	Immunogen	Recombinant Human APOD protien	Recombinant Human APOD protien	
Information	Swissprot	P05090		
Product details	Reactivity	Human	Human	
	Host	Rabbit	Rabbit	
	Conjugation	Unconjugated	Biotin	
	Concentration	0.5 mg/mL	/	
	Buffer	PBS with 0.04% Proclin 300; 50%	PBS with 0.04% Proclin 300; 1%	
		glycerol; pH 7.5	protective protein; 50% glycerol; pH	
			7.5	
	Purify	Antigen Affinity	Antigen Affinity	
	Specificity	Detects Human APOD in ELISAs.		

For Research Use Only

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017 Web: www.elabscience.com Email: techsupport@elabscience.com



A Reliable Research Partner in Life Science and Medicine

Applications

Human APOD Sandwich ELISA Assay

	Recommended	Reagent	Images
	Concentration/Dilution		
ELISA	0.5-4 μg/mL	Human ApoD Capture	
Capture		Antibody	10
			usity 13
ELISA	1:1000-1:10000	Human ApoD Detection	Optical Density
Detection	111000 1110000	Antibody (Biotin)	0.1
			1 10 100 100 1000 Human ApoD Concentration (ng/mL)

Note: This standard curve is only for demonstration purposes. A standard curve should be generated for each assay!

Background

Apolipoprotein D (ApoD) is a member of the lipocalin superfamily of ligand transporters; and has been implicated in the transport of small hydrophobic molecules. ApoD is also a component of plasma high-density lipoproteins (HDL) . Alteration of ApoD expression has been linked to multiple neurological disorders; including Alzheimer's disease.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017 Web: www.elabscience.com Email: techsupport@elabscience.com