

The professional cell culture empowers a healthier world

rocellsystem

DMEM (Glucose free)

Cat. No.: PM150270

Size: 500mL

General Information

Product Form Liquid Concentration 1 × 7.2-7.4 pΗ Negative (R **D-Glucose** Negative **HEPES** L-Glutamine 4 mM NaHCO₃ 3700 mg/L

Phenol red 15 mg/L Sodium pyruvate 1 mM

2-8°C, Shading Light Storage

Room Temperature **Shipping**

Expiration date 24 months

Background

DMEM (Dulbecco's Modified Eagle Medium) was developed on the basis of MEM medium. Compared with MEM medium, the content of amino acid increased by 2 times, the content of vitamin increased by 4 times, and the content of non-essential amino acid, trace iron ion and sodium pyruvate were increased by 4 times.

The glucose content of DMEM medium was originally designed as 1000 mg/L (low Glucose type), and then developed into 4500 mg/L (high Glucose type), which has been widely used in cell culture.

DMEM medium without glucose can be added to the concentration of glucose at will according to the research needs, which is convenient and quick.

DMEM (glucose free) contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors, so the product should be used with serum or serum-free additives. Procell system

Notes

- This product is for research use only.
- This product is sterilized by 0.1 µm filtration.
- It is necessary to pay attention to the aseptic operation and avoid the contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Web: www.procellsystem.com Email: techsupport@procellsystem.com