

The professional cell culture empowers a healthier world

ocell system

all system

## DMEM (Low glucose), without phenol red

Cat. No.: PM150251

Size: 500mL

## **General Information**

Product Form

Concentration

pH

7.2-7.4

D-Glucose

HEPES

L-Glutamine

Liquid

1 ×

7.2-7.4

1000 mg/L

Negative

4 mM

NaHCO₃3700 mg/LPhenol redNegativeSodium pyruvate1 mM

**Storage** 2-8°C, Shading Light

Shipping Room Temperature

**Expiration date** 24 months

## **Background**

DMEM (Dulbecco's Modified Eagle Medium) was developed on the basis of MEM medium. Compared with MEM medium, the content of amino acid increased by 2 times, the content of vitamin increased by 4 times, and the content of non-essential amino acid, trace iron ion and sodium pyruvate were increased by 4 times.

The glucose content of DMEM medium was originally designed as 1000 mg/L (low Glucose type), and then developed into 4500 mg/L (high Glucose type), which has been widely used in cell culture.

DMEM (Low glucose) contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors, so the product should be used with serum or serum-free additives.

## **Notes**

1. This product is for research use only.

Web: www.procellsystem.com

- This product is sterilized by 0.1 µm filtration.
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Email: techsupport@procellsystem.com