

Ham's F-12, with P/S

Cat. No.: PM150810A

Size: 500mL

General Information

Product Form Liquid Concentration 1 × 7.2 - 7.4pН **D-Glucose** 1802 mg/L **Sodium Bicarbonate (NaHCO3)** 1176 mg/L L-Glutamine 1 mM **HEPES** Negative Phenol Red 1.2 mg/L **Sodium Pyruvate** 1 mMStreptomycin Sulfate 100 mg/L Penicillin G 100 kU/L

Pricella by Elabscience 2-8°C, Shading Light Storage Room Temperature **Shipping**

Expiration date 24 months

Background

Ham 's F-12 Nutrient Mixture was designed by Ham in 1969 based on Ham 's F-10 nutrient mixture, which was originally used for serum-free culture of CHO cells. Ham 's F-12 is often used as the basic culture medium for serum free culture. At low serum levels, it is especially suitable for single cell culture and clonal culture. After adding serum, it is also widely used in the culture of cancer cells and primary cells, such as rat hepatocytes, rat prostate epithelial cells, chondrocytes, rat myoblasts, and chicken embryonic cells. In addition, the DMEM/F12 medium obtained by mixing Ham 's F-12 in equal volume with DMEM is more nutritious and widely used. icella

Notes

- 1. This product is for research use only.
- 2. This product is sterilized by 0.1 µm filtration.
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Web: www.pri-cella.com Email: techsupport@pri-cella.com