

Normal Goat Blocking Buffer (Ready-to-Use)



Cat. No: E-IR-R110

Size: 5 mL / 10 mL / 50 mL

Product Content

Cat.	Products	5 mL	10 mL	50 mL	Storage
E-IR-R110	Normal Goat Blocking Buffer (Ready-to-Use)	5 mL	10 mL	50 mL	2~8 °C

Introduction

This product is normal goat Blocking Buffer, added stabilizer, As a blocking buffer, It can effectively reduce the nonspecific staining of IHC/IF experiment, so as to reduce the background.

Instructions

This product is ready to use goat Blocking Buffer and can be used directly.

When using as control IgG, the dosage should be the same as that of other IgG. Too much of the control IgG may cause the background to be too high.

When using as an immunohistochemical blocking buffer, please directly drop it onto the tissue slices for blocking.

Storage

Store at 2~8 °C for 12 months.

Cautions

1. This product couldn't be used for clinical diagnosis or treatment, food or medicine, and can't be stored in residence.
2. For your safety and health, please wear the lab coat and disposable gloves before the experiments.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Rev. V1.3