

The professional cell culture empowers a healthier world

cellsystem

MFMa

Cat. No.: PM150421

Size: 500mL

General Information

Product Form Liquid Concentration 1 × Hq 7.2 - 7.41000 mg/L **D-Glucose** sell sy Negative **HEPES** 2 mM L-Glutamine NaHCO₃ 2200 mg/L

Phenol red 10 mg/L Sodium pyruvate 1 mM **Nucleoside** Positive **Balanced salt** Earle's salt

rocell system Storage 2-8°C, Shading Light **Shipping** Room Temperature

Expiration date 24 months

Background

MEMα was developed on the basis of MEM medium. Compared with MEM medium, NEAA, sodium pyruvate, zinc sulfate, VB12, biotin and ascorbic acid were added. It is widely used in suspension and adherent cell culture of various mammals. Mem a medium without nucleoside and deoxynucleoside is often used as medium for DG44 and other DHFR- cells.

This product contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors. Therefore, the product should be used with serum or serum-free additives.

Notes

- This product is for research use only.
- This product is sterilized by 0.1 µm filtration.
- It is necessary to pay attention to the aseptic operation and avoid the contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017 Web: www.procellsystem.com

Email: techsupport@procellsystem.com