

## Recombinant Human MFAP4 Protein (His Tag)

**Catalog Number:** PKSH032754

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


### Description

|                      |   |
|----------------------|---|
| <b>Species</b>       | Human   |
| <b>Source</b>        | HEK293 Cells-derived Human MFAP4 protein Val22-Ala255, with an C-terminal His |
| <b>Calculated MW</b> | 27.5 kDa  |
| <b>Observed MW</b>   | 37 kDa  |
| <b>Accession</b>     | AAH62415.1  |
| <b>Bio-activity</b>  | Not validated for activity  |

### Properties

|                       |  |
|-----------------------|--|
| <b>Purity</b>         | > 85 % as determined by reducing SDS-PAGE.   |
| <b>Endotoxin</b>      | < 1.0 EU per µg of the protein as determined by the LAL method.  |
| <b>Storage</b>        | Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months. |
| <b>Shipping</b>       | This product is provided as lyophilized powder which is shipped with ice packs.  |
| <b>Formulation</b>    | Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization.   |
| <b>Reconstitution</b> | Please refer to the specific buffer information in the printed manual.   |

### Data



> 85 % as determined by reducing SDS-PAGE.

### Background

Microfibril-associated glycoprotein 4(MFAP4) is a secreted protein and contains 1 fibrinogen C-terminal domain, similarity to a bovine microfibril-associated protein. The protein has binding specificities for both collagen and carbohydrate. It is thought to be an extracellular matrix protein which is involved in calcium-dependent cell adhesion or intercellular interactions. The gene is located within the Smith-Magenis syndrome region.

### For Research Use Only

Toll-free: 1-888-852-8623  
Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086  
Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017