he professional cell culture enpowers a healthier world

DMEM (Low glucose), without sodium pyruvate, phenol red

Cat. No: PM150253

Size: 500mL

General Information

Product Form

Concentration

pH

7.2-7.4

D-Glucose

HEPES

L-Glutamine

NaHCO3

Liquid

1×

1000mg/L

Negative

4mM

3700mg/L

Phenol red Negative
Sodium pyruvate negative

Storage 2-8℃, Shading Light

Shipping RT

Expiration date 24 months



Background

DMEM (Dulbecco's Modified Eagle Medium) was developed on the basis of MEM medium. Compared with MEM medium, the content of amino acid increased by 2 times, the content of vitamin increased by 4 times, and the content of non-essential amino acid, trace iron ion and sodium pyruvate were increased by 4 times. The glucose content of DMEM medium was originally designed as 1000 mg/L (low Glucose type), and then developed into 4500mg/L (high Glucoser type), which has been widely used in cell culture.

DMEM (Low glucose) contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors, so the product should be used with serum or serum-free additives.

Notes

Web: www.pri-cella.com

- 1. This product is for research use only.;
- 2. This product is sterilized by 0.1µm filtration.;
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Email: techsupport@pri-cella.com