

Vero Serum-Free Medium

Cat. No.: SF4001

Size: 500mL

Description

Vero Serum-Free Medium is a chemically defined medium designed and developed specifically for Vero cells. It is suitable for Vero cell adherent culture. Vero Serum-Free Medium contains L-Glutamine, which generally does not require additional addition.

General Information

Product Form	Liquid
Sterile Filtered	Yes
Specification	500mL
pH	7.2-7.6
Storage	2-8°C, Shading Light
Shipping	Ice bag
Expiration date	12 months

Use Instructions

This product can be used directly for serum-free cultured Vero cells. If the culture is transferred from serum-containing medium to a serum-free medium, it is necessary to domesticate Vero cells using a serum-free medium. The taming gradient is 2%-1%-0.5%-0% (the proportion of serum added to Vero Serum-Free Medium), and serum-free cell domestication can be completed within 1-2 weeks.

Notes

1. This product is only used for scientific research or further research, not for diagnosis and treatment.
2. This product has been filtered and sterilized. Pay attention to aseptic operation to avoid contamination.
3. To maintain the best use of this product, do not perform freeze-thaw treatment.