

The professional cell culture empowers a healthier world

HAT Solution, 50 ×

Cat. No.: PB180631

Size: 10mL

General Information

Product Form Liquid Concentration 50× **Hypoxanthine** 5 mM 0.02 mM **Aminopterin** Cell 5 y 0.8 mM **Thymidine** Dissolvent 0.85% NaCl Storage -5~-20°C **Shipping** Ice bag **Expiration date** 12 months

Background

Procell system HAT solution is a mixture of hypoxanthine, aminopterin and thymidine. It is often used for the screening and culture of hybridoma cells.

Notes

- This product is for research use only.
- This product is sterilized by 0.1 µm filtration.
- It is necessary to pay attention to the aseptic operation and avoid the pollution.
- This product is 50 × concentrated solution and dilute it to to the desired concentration before use.





Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Email: techsupport@procellsystem.com Web: www.procellsystem.com