

The professional cell culture empowers a healthier world

ocell system

ellsystem

DMEM (Low glucose)

Cat. No. : PM150220 Size : 500mL / 1000mL

General Information

Product Form Liquid

Concentration 1 ×

pH 7.2-7.4

D-Glucose 1000 mg/L

Sodium Bicarbonate (NaHCO₃) 3700 mg/L

L-Glutamine 4 mM
HEPES Negative
Phenol Red 15 mg/L
Sodium Pyruvate 1 mM

Storage 2-8°C, Shading Light

Shipping Room Temperature

Expiration date 24 months

Background

DMEM (Dulbecco's Modified Eagle Medium) was developed on the basis of MEM medium. Compared with MEM medium, the content of amino acid increased by 2 times, the content of vitamin increased by 4 times, and the content of non-essential amino acid, trace iron ion and sodium pyruvate were increased by 4 times.

The glucose content of DMEM medium was originally designed as 1000 mg/L (low Glucose type), and then developed into 4500 mg/L (high Glucoser type), which has been widely used in cell culture. DMEM (Low glucose) contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors, so the product should be used with serum or serum-free additives.

Notes

- 1. This product is for research use only.
- This product is sterilized by 0.1 μm filtration.
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Web: www.procellsystem.com Email: techsupport@procellsystem.com