



A Reliable Research Partner in Life Science and Medicine

Rat HPX Antibody Pair Set

Catalog No. E-KAB-0392 Applications ELISA

Synonyms HPX, HX, hemopexin, Beta-1B-glycoprotein

Kit components & Storage

Title	Specifications	Storage
Rat HPX Capture Antibody	1 vial, 100 μ g	Store at -20°C for one year.
		Avoid freeze / thaw cycles.
Rat HPX Detection Antibody (Biotin)	1 vial, 50 μL	Store at -20°C for one year.
		Avoid freeze / thaw cycles.

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Items		Characteristic (E-KAB-0392)		
		Rat HPX Capture Antibody	Rat HPX Detection Antibody (Biotin)	
Immunogen	Immunogen	Recombinant Rat HPX protein	Recombinant Rat HPX protein	
Information	Swissprot	P20059		
Product details	Reactivity	Rat	Rat	
	Host	Rabbit	Rabbit	
	Conjugation	Unconjugated	Biotin	
	Concentration	0.5mg/mL	/	
	Buffer	PBS with 0.04% Proclin 300, 50%	PBS with 0.04% Proclin 300, 1%	
		glycerol, pH 7.4	protective protein, 50% glycerol, pH	
			7.4	
	Purify	Protein A	Protein A & Antigen Affinity	
Specificity		Detects Rat HPX in ELISAs.		

For Research Use Only

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017 Web: www.elabscience.com Email: techsupport@elabscience.com



Elabscience®

A Reliable Research Partner in Life Science and Medicine

Applications

Rat HPX Sandwich ELISA Assay:

	Recommended	Reagent	Images
	Concentration/Dilution		
ELISA	0.5-4µg/mL	Rat HPX Capture Antibody	
Capture			Nigs 1
ELISA Detection	1:1000-1:10000	Rat HPX Detection Antibody (Biotin)	0.01 100 1000 10000 Rat HPX concentration(ng/mL)

Note: This standard curve is only for demonstration purposes. A standard curve should be generated for each assay!

Background

Hemopexin (HPX) is the plasma protein responsible for scavenging heme, thus preventing heme-mediated oxidative stress and heme-bound iron loss. In addition, hemopexin blocks heme activation of immune receptors and vascular inflammatory processes. It is mainly expressed in liver, the synthesis of which is induced after inflammation. Alterations of plasma hemopexin level have been linked to disorders like atherosclerosis and inflammatory diseases.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017 Web: www.elabscience.com Email: techsupport@elabscience.com