

The professional cell culture empowers a healthier world

IMDM

Cat. No.: PM150510

Size: 500mL

General Information

Product Form Liquid Concentration 1 × pΗ 7.2 - 7.44500 mg/L **D-Glucose** cell sy 25 mM **HEPES** L-Glutamine

NaHCO₃ 3024 mg/L Phenol red 15 mg/L Sodium pyruvate 1 mM

rocellsystem 2-8°C, Shading Light **Storage**

Room Temperature **Shipping**

Expiration date 24 months

Background

IMDM (Iscove's Modified Dulbecco Medium)is an improved DMEM medium, which is used for the culture of erythrocyte progenitor cells and macrophages. The IMDM medium was supplemented with selenium, HEPES, sodium pyruvate and additional amino acids and vitamins on the basis of the DMEM medium, and potassium nitrate instead of iron nitrate, which is suitable for the rapid proliferation, high-density cell culture. IMDM medium can not only culture cells with special nutritional requirements (such as mouse B lymphocytes, stimulated B cells stimulated by LPS, bone marrow hematopoietic cells, T lymphocytes and various hybridoma cells), It can also be used as a base solution for some unique serum-free medium.

This product contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors. Therefore, the product should be used with serum or serum-free additives. Procell system

Notes

- This product is for research use only.
- This product is sterilized by 0.1 µm filtration. 2.
- It is necessary to pay attention to the aseptic operation and avoid the contamination.

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Web: www.procellsystem.com Email: techsupport@procellsystem.com