## IL-4/BSF-1, Human, Recombinant

Cat. No. : PCK031

产品信息	
别名	BCGF;BCDF;B-cell Stimulating Factor (BSF-1)
物种	Human rocell system
表达宿主、日本主	E.coli
序列信息	MHKCDITLQEIIKTLNSLTEQKTLCTELTVTDIFAASKNTTEKETFCRAATVLRQFYSHHEKD
E	$\label{eq:transf} TRCLGATAQQFHRHKQLIRFLKRLDRNLWGLAGLNSCPVKEANQSTLENFLERLKTIMREKY$
	SKCSS with polyhistidine tag at the C-terminus.
检索号	P05112.1
标签	His-tag at the C-terminus
分子量	His-tag at the C-terminus 15.9 kDa 12 months Macana ha it a hilita to in the TTP 1 HIL - Liferation The TD 50 for this offset is 6.02 m/mL
有效期	12 months ® Proceedings
生物活性	Measure by its ability to induce 1F-1 cens promeration. The ED50 for this effect is < 0.2 ng/mL.
	The specific activity of recombinant human IL-4 is approximately $> 2.8 \times 10^{7}$ IU/mg.
产品特性	
内毒素	< 0.1 EU per 1 µg of the protein by the LAL method.
保存	Lyophilized protein should be stored at -5~-20°C for 1 year. Upon reconstitution, store at 2-8°C
	for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1%
	BSA, 10% FBS, 5% HSA or 5% trehalose solution), protein aliquots should be stored at -5~-20
\= + <u>A</u>	°C or -80°C for 3-6 months.
运输	Ambient temperature or ice pack.
制剂	The protein was lyophilized from a 0.2 $\mu$ m filtered solution containing 1 × PBS, pH 8.0.
复溶	It is recommended to reconstitute the lyophilized protein in sterile water to a concentration not
E CE	less than 100 µg/mL. Do Not Vortex! Vigorous shaking may impair the biological activity of the
背景介绍	protein.

The interleukin 4 (IL4, IL-4) is a cytokine that induces differentiation of naive helper T cells (Th0 cells) to Th2 cells. Upon activation by IL-4, Th2 cells subsequently produce additional IL-4 in a positive feedback loop. The cell that initially produces IL-4, thus inducing Th2 differentiation, has not been identified, but recent studies suggest that basophils may be the effector cell. It is closely related and has functions similar to Interleukin 13.



## 普诺赛<sup>®</sup> Procell system



普诺赛<sup>®</sup>\Procellsystem



普诺赛<sup>®</sup>\Procellsystem



Rev. V1.2