

A Reliable Research Partner in Life Science and Medicine

Serum Albumin Monoclonal Antibody(Detector)

catalog number: AN001410P

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity Human

Immunogen Human ALB Native Protein

Host Mouse

Isotype Mouse IgG2b

Clone 9C5


Purification Protein A/G Purification

Buffer Phosphate buffered solution, pH 7.2, containing 0.05% proclin 300.

Applications Recommended Dilution

ELISA Detector $0.1-0.4 \mu g/mL$

Data



Sandwich ELISA-Human ALB Native Protein standard curve. Background subtracted standard curve using Serum Albumin antibody(AN001400P)(Capture), Serum Albumin antibody(AN001410P)(Detector) in sandwich ELISA. The reference range value for Human ALB Native Protein is 0.781-50 ng/mL.

Preparation & Storage

Storage Storage Store at 4°C valid for 12 months or -20°C valid for long term storage, avoid freeze /

thaw cycles.

Shipping The product is shipped with ice pack, upon receipt, store it immediately at the

temperature recommended.

Background

For Research Use Only

Fax: 1-832-243-6017

Elabscience Bionovation Inc.



A Reliable Research Partner in Life Science and Medicine

Binds water,Ca2+,Na+,K+,fatty acids,hormones,bilirubin and drugs (Probable). Its main function is the regulation of the colloidal osmotic pressure of blood (Probable). Major zinc transporter in plasma,typically binds about 80% of all plasma zinc.Major calcium and magnesium transporter in plasma,binds approximately 45% of circulating calcium and magnesium in plasma (By similarity).

Potentially has more than two calcium-binding sites and might additionally bind calcium in a non-specific manner (By similarity). The shared binding site between zinc and calcium at residue Asp-273 suggests a crosstalk between zinc and calcium transport in the blood (By similarity). The rank order of affinity is zinc > calcium > magnesium (By similarity). Binds to the bacterial siderophore enterobactin and inhibits enterobactin-mediated iron uptake of E.coli from ferric transferrin, and may thereby limit the utilization of iron and growth of enteric bacteria such as E.coli. Does not prevent iron uptake by the bacterial siderophore aerobactin

Toll-free: 1-888-852-8623 Web:www.elabscience.com

Tel: 1-832-243-6086 Email:techsupport@elabscience.com