

Recombinant Human SUS D4/Sushi domain-containing protein 4 Protein (Fc Tag)



Catalog Number:PKSH030623

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	PRO222
Species	Human
Expression Host	HEK293 Cells
Sequence	Met 1-Phe290
Accession	Q5VX71-3
Calculated Molecular Weight	53.8 kDa
Observed molecular weight	67 kDa
Tag	C-mFc

Properties

Purity	> 85 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4 Normally 5 % - 8 % trehalose, mannitol and 0.01 % Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 85 % as determined by reducing SDS-PAGE.

Background

SUSD4, also known as sushi domain-containing protein 4, is a hypothetical cell surface protein whose tissue distribution and function are completely unknown. SUSD4 is detectable in murine brains, eyes, spinal cords, and testis but not other tissues. In brains, SUSD4 is highly expressed in the white matter on oligodendrocytes/axons, and in eyes, it is exclusively expressed on the photoreceptor outer segments. In vitro complement assays, SUSD4 augments the alternative but not the classical pathway of complement activation at the C3 convertase step. SUSD4 deficiency may cause autism or Fryns syndrome, both of which are genetic diseases with severe abnormal neurological development and/or functions.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017