



A Reliable Research Partner in Life Science and Medicine

FITC Armenian Hamster IgG Isotype Control[PIP]

Catalog Number: E-AB-F09852C

Note: Centrifuge before opening to ensure complete recovery of vial contents.

| _ | | | | - | | 40 | | |
|----|-----|---|---|---|---|-----|---------------------|---|
| | e: | e | ^ | 7 | | TI. | $\boldsymbol{\cap}$ | n |
| ш, | 4-1 | - | • | | w | ш | u | |

Host Armenian Hamster Isotype Armenian Hamster IgG

Clone No. PIP
Conjugation FITC

Conjugation Information FITC is designed to be excited by the Blue laser (488 nm) and detected using an optical

filter centered near 530 nm (e.g., a 525/40 nm bandpass filter).

Storage Buffer Phosphate buffered solution, pH 7.2, containing 0.09% sodium azide and 1% BSA.

Applications Recommended usage

FCM Each lot of this antibody is quality control tested by flow cytometric analysis. The amount

Web: www.elabscience.cn

of the reagent is suggested to be used 5 μ L of antibody per test (million cells in 100 μ L staining volume or per 100 μ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for

individual use.

Preparation & Storage

Storage Keep as concentrated solution.

This product can be stored at 2-8°C for 12 months. Please protected from prolonged

exposure to light and do not freeze.

Shipping Ice bag