ocellsystem





MEM/F12, without phenol red

Cat. No.: PM151226

Size: 500mL

General Information

Product Form Liquid 1 × Concentration pН 7.2-7.4 1401 mg/L **D-Glucose** 1688 mg/L Sodium Bicarbonate (NaHCO₃) L-Glutamine 1.5 mM **Phenol Red** Negative **Balanced salt** Earle's salt **Sodium Pyruvate** $0.5 \, \text{mM}$

2-8°C, shading light Storage

Room Temperature Shipping

Expiration date 24 months

Background

MEM/F12 medium is based on MEM medium with more abundant nutrients in F12 medium and contains a variety of trace elements. It is widely used in the cultivation of a variety of mammalian cells. At the same time, MEM/F12 medium is often used as the basis for the development of serum-free medium, and it is also suitable for the cultivation of mammalian cells with low serum content. This product contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors. Therefore, the product should be used with serum or serum-free additives.

Notes

Toll-free: 1-888-852-8623

- This product is for research use only.
- 2. This product is sterilized by 0.1 µm filtration.
- It is necessary to pay attention to the aseptic operation and avoid the contamination. 3. procellsy

Tel: 1-832-243-6086 Fax: 1-832-243-6017 Web: www.procellsystem.com

Email: techsupport@procellsystem.com