

## MEMα, with L-alanyl-L-glutamine, without nucleoside

Liquid

Cat. No: PM150426

Size: 500mL

**Product Form** 

## **General Information**

Concentration 7.2-7.4 Ηα 1000mg/L **D-Glucose** 

**HEPES** Negative L-Alanyl-L-Glutamine 2mM NaHCO3 2200mg/L Phenol red 10mg/L Sodium pyruvate 1 mM **Balanced salt** Earle's salt **Nucleoside** Negative

Storage 2-8℃, Shading Light

**Shipping** RT

**Expiration date** 24 months

## **Background**

MEM α was developed on the basis of MEM medium. Compared with MEM medium, NEAA, sodium pyruvate, zinc sulfate, VB12, biotin and ascorbic acid were added. It is widely used in suspension and adherent cell culture of various mammals. Mem  $\alpha$  medium without nucleoside and deoxynucleoside is often used as medium for DG44 and other DHFR- cells.

This product contains many kinds of amino acids, vitamins, inorganic salts and other ingredients for cell culture, but does not contain protein, lipids or any growth factors. Therefore, the product should be used with serum or serum-free additives.

## Notesby Elab

- 1. This product is for research use only.;
- 2. This product is sterilized by 0.1µm filtration.;
- 3. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture. by Elabscie

Email: techsupport@pri-cella.com

Pricella®
by Elabscience