

A Reliable Research Partner in Life Science and Medicine

Recombinant Histone H3 (Mono Methyl Lys122) Monoclonal Antibody

catalog number: AN302115L

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity Human; Rat; Mouse

Immunogen Monomethylated human histone H3 (Lys122) peptide

 Host
 Rabbit

 Isotype
 IgG, κ

 Clone
 A839

Purification Protein A purified

Buffer PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protein protectant.

Applications Recommended Dilution

WB 1:2000-1:10000

Data



Western Blot with Histone H3 (Mono Methyl Lys122)

Monoclonal Antibody at dilution of 1:10000. Lane 1: HeLa.

Lane 2: 4T1, Lane 3: BRL, Lane 4: Mouse liver, Lane 5: Rat

brain, Lane 6: Recombinant Histone H3 (20ng)

Observed-MW:15 kDa Calculated-MW:15 kDa

Preparation & Storage

Storage Storage Store at -20°C Valid for 12 months. Avoid freeze / thaw cycles.

Shipping Ice bag

Background

Histone post-translational modifications (PTMs) are key mechanisms of epigenetics that modulate chromatin structures, termed as "histone code". The PTMs on histone including acetylation, methylation, phosphorylation and novel acylations directly affect the accessibility of chromatin to transcription factors and other epigenetic regulators, altering genome stability, gene transcription, etc. Histone methylation occurs primarily at lysine and arginine residues on the amino terminal of core histones. Methylation of histones can either increase or decrease transcription of genes, depending on which amino acids (Lys or Arg) in the histones are methylated and how many methyl groups are attached (mono-, di-, Trimethylation on Lys, mono-di-symmetric/asymmetric methylation on Arg). Mostly, lysine methylation occurs primarily on histone H3 Lys4, 9, 27, 36, 79 and H4 Lys20, while Arginine methylation occurs primarily on histone H3 Arg2, 8, 17, 26 and H4 Arg3. histone methyltransferases (HMTs) and histone demethylases (HDMs) are major regulating factors.

For Research Use Only

 Toll-free: 1-888-852-8623
 Tel: 1-832-243-6086
 Fax: 1-832-243-6017

 Web: www.elabscience.com
 Email: techsupport@elabscience.com
 Rev. V1.0